


Arc Tipping Lever



Not able to retrofit

Item Number	Part Number	Description	Quantity
1 - 7	120814	Tipping Lever Assembly Arc	1
3, 5 - 7	120815	Tipping Lever Mounting Hardware Arc	1
3	110269	Mount Extension, Vent Tray / Tipping Lever	2
4	100804	Tube End Plug 1"	2
5	101810	M6 x 35 SHCS BLK ZC	2

Item Number	Part Number	Description	Quantity
6	101847	M6 x 40 SHCS Blk Zn	2
7	100658	M6 Nylock Nut Blk Zn	4