


Ride Designs® Corbac™ Adjustable Contour Back Support



The unique spacer fabric material is used in both the Corbac and Forward Cushion for improved comfort and breathability.



▲ First, dial in the exact posterior contour with the Corbac's adjustment straps.

Then, refine the fit even further with the thoraco-lumbar-pelvic (TLP) supports. TLP supports provide precise adjustment of posterior and lateral contours. ▶

Lightweight

At two pounds, the Corbac™ weighs little more than standard back upholstery.

Precision fit

The Corbac is easily adjusted to help improve postural control and alignment. It does not decrease the wheelchair's depth or width. The sitter experiences fitting in, not sitting on, the wheelchair.



Reduce heat and moisture buildup

The spacer fabric overlay promotes air circulation and wicks moisture from the back surface.

Versatile

No screws, nuts, or tricky hardware, just a quick and easy mounting system that takes seconds to assemble onto the wheelchair.

A cinch to clean

The overlay detaches easily and can be machine washed.