

LIBERTY FT

US Pat# 10,828,213



LIBERTY IS INDEPENDENCE

Finally, one wheelchair that offers true portability, enables independent propulsion and provides the benefits of tilt-in-space!



PROPULSION

Liberty is designed for independent mobility, especially foot propulsion. The front seat height remains the same regardless of the tilt angle. Simply choose the tilt angle that provides the best pelvic stability without compromising mobility.

REPOSITIONING

Everyone needs relief from sitting in a static position. Whether your primary concern is pressure relief, facilitating proper posture or just providing a break from gravity, Liberty makes it painless. Featuring an easy to adjust gas spring design, adjusting the system for effortless repositioning is simple.



PORTABILITY

Liberty FT is the only adult tilt-in-space that can fold into a portable package ideal for transportation. Want an even smaller package? Choose the optional folding backrest. With a transport weight under 26 pounds, it's the lightest adult tilt-in-space available. Liberty FT gives you the independence to get out and go!



OPTION TOUR



Liberty is one of the most versatile solutions on the market providing the options needed to meet the end users changing needs. In most cases, the modular design allows for easy change of dimensional characteristics making it the ideal choice for short discharge notification in stroke rehabilitation.

ACCESSORIES

O₂ Holder
IV Pole
Transit Wheels

ARMRESTS

Height Adjustable
Angle Adjustable Locking
Extendable Flip-Up

SEAT & BACK

Growing Seat Pan
Standard Upholstery
Tension Adjustable
Axiom Back

BACK CANES

Stroller Handle
Height Adjustable
Height Adjustable with 8° Bend
Angle Adjustable (Folding)
Fixed (Non-Folding)

CASTERS

Widths: 1" | 1.5" | 2"
Heights: 4" | 5" | 6" | 7" | 8"

FOOTPLATES

Composite
Composite Angle Adjustable
Aluminum Angle Adjustable
Aluminum Locking Multi-Angle Adjustable
One Piece Flip-Up Angle Adjustable

HANGERS

Front Mount
Extension Tube
Pro ELR
Pediatric Pro ELR

REAR WHEELS

12" | 16" | 20" | 22" | 24"

WHEEL LOCKS

Push to Lock
Pull to Lock
Hemi Wheel Lock

SPECIFICATIONS

Tilt Range:	0° - 20°	0° - 20°
Seat Width:	14 - 22"	36 - 56 cm
Seat Depth:	14 - 20"	36 - 51 cm
Seat Heights:	13 - 20"	33 - 51 cm
Liberty FT Transport Weight *:	25.8 lbs.	11.7 kg
Weight Capacity:	250 lbs.	113 kg
WC-19 Approved **:	Yes	Yes

* 16x16 fixed back with armrests, footrests and rear wheels removed.

** Liberty is WC-19 approved with the folding backrest only.



Stevens Point, WI | 800-981-1540 | www.kimobility.com

DCN0314.4