

## CR45 Wheels - 16" Mag



*If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. if changing wheel size and Attendant Foot Lock, spacer 112503 is required which can be found on the CR45 Axle Plate page*

Item Number	Part Number	Description	Quantity
1a	001159	16" Mag Wheel w/ Pneumatic Assembly Tilt	1
1b	001160	16" Mag Wheel w/ Pneumatic w/ Airless Insert Assembly Tilt	1
2	001161	16" Mag Wheel w/ Full Poly Assembly Tilt	1
3	003453	16" Mag Wheel w/ Low Poly Assembly Tilt	1