

Ethos Handrims - Natural Fit



For more information, see Handrim Hardware Chart. Quantity in Handrim Hardware Kits are for one wheel. <p>Tab mounted handrims must be installed prior to the tire being installed.

Item Number	Part Number	Description	Quantity	Retail
1a	200538	Handrim Natural Fit Standard 6 Rivet Nut 22" <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i>	1	\$213.41
1a	107168	Handrim Natural Fit Standard 6 Rivet Nut 22" w/ Micro Grip <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i>	1	\$236.43
1b	100793	Handrim Natural Fit Standard 6 Rivet Nut 24" <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i>	1	\$213.41
1b	107169	Handrim Natural Fit Standard 6 Rivet Nut 24" w/ Micro Grip <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i>	1	\$251.11
1c	200539	Handrim Natural Fit Standard 6 Rivet Nut 25" <i>Used on 25" Ki Spoke and 25" MAXX Spoke</i>	1	\$213.41
1c	107170	Handrim Natural Fit Standard 6 Rivet Nut 25" w/ Micro Grip <i>Used on 25" Ki Spoke and 25" MAXX Spoke</i>	1	\$226.40
1d	100907	Handrim Natural Fit Standard 6 Rivet Nut 26" <i>Used on 26" Ki Spoke and 26" MAXX Spoke</i>	1	\$213.41
1d	107171	Handrim Natural Fit Standard 6 Rivet Nut 26" w/ Micro Grip <i>Used on 26" Ki Spoke and 26" MAXX Spoke</i>	1	\$255.11
1e	100768	Handrim Natural Fit Standard 5 Rivet Nut 24" <i>Used on 24" 5-Spoke Mag</i>	1	\$213.41
1e	107172	Handrim Natural Fit Standard 5 Rivet Nut 24" w/ Micro Grip <i>Used on 24" 5-Spoke Mag</i>	1	\$239.15

Item Number	Part Number	Description	Quantity	Retail
1f	100889	Handrim Natural Fit Standard 6 Tab 22" <i>Used on 22" Spinergy</i>	1	\$213.41
1f	107173	Handrim Natural Fit Standard 6 Tab 22" w/ Micro Grip <i>Used on 22" Spinergy</i>	1	\$236.43
1g	100830	Handrim Natural Fit Standard 6 Tab 24" <i>Used on 24" Spinergy/Superlight/Topolino</i>	1	\$213.41
1g	107174	Handrim Natural Fit Standard 6 Tab 24" w/ Micro Grip <i>Used on 24" Spinergy /Superlight/Topolino</i>	1	\$239.15
1h	101464	Handrim Natural Fit Standard 6 Tab 25" <i>Used on 25" Spinergy/Superlight/Topolino</i>	1	\$213.41
1h	107175	Handrim Natural Fit Standard 6 Tab 25" w/ Micro Grip <i>Used on 25" Spinergy/Superlight/Topolino</i>	1	\$240.24
1i	200200	Handrim Natural Fit Standard 6 Tab 26" <i>Used on 26" Spinergy</i>	1	\$234.83
1i	107176	Handrim Natural Fit Standard 6 Tab 26" w/ Micro Grip <i>Used on 26" Spinergy</i>	1	\$263.36
2a	200201	Handrim Natural Fit LT 6 Rivet Nut 22" <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i>	1	\$234.83
2a	107177	Handrim Natural Fit LT 6 Rivet Nut 22" w/ Micro Grip <i>Used on 22" Ki Spoke and 22" MAXX Spoke</i>	1	\$256.83
2b	200202	Handrim Natural Fit LT 6 Rivet Nut 24" <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i>	1	\$234.83
2b	107178	Handrim Natural Fit LT 6 Rivet Nut 24" w/ Micro Grip <i>Used on 24" Ki Spoke and 24" MAXX Spoke</i>	1	\$260.64
2c	200540	Handrim Natural Fit LT 6 Rivet Nut 25" <i>Used on 25" Ki Spoke and 25" MAXX Spoke</i>	1	\$234.83
2c	107179	Handrim Natural Fit LT 6 Rivet Nut 25" w/ Micro Grip <i>Used on 25" Ki Spoke and 25" MAXX Spoke</i>	1	\$273.67
2d	101454	Handrim Natural Fit LT 6 Rivet Nut 26" <i>Used on 26" Ki Spoke and 26" MAXX Spoke</i>	1	\$234.83
2d	107180	Handrim Natural Fit LT 6 Rivet Nut 26" w/ Micro Grip <i>Used on 26" Ki Spoke and 26" MAXX Spoke</i>	1	\$263.36
2e	100769	Handrim Natural Fit LT 5 Rivet Nut 24" <i>Used on 24" 5-Spoke Mag</i>	1	\$234.83
2e	107181	Handrim Natural Fit LT 5 Rivet Nut 24" w/ Micro Grip <i>Used on 24" 5-Spoke Mag</i>	1	\$259.55

Item Number	Part Number	Description	Quantity	Retail
2f	100888	Handrim Natural Fit LT 6 Tab 22" <i>Used on 22" Spinergy</i>	1	\$234.83
2f	107182	Handrim Natural Fit LT 6 Tab 22" w/ Micro Grip <i>Used on 22" Spinergy</i>	1	\$269.67
2g	100828	Handrim Natural Fit LT 6 Tab 24" <i>Used on 24" Spinergy/Superlight/Topolino</i>	1	\$234.83
2g	107183	Handrim Natural Fit LT 6 Tab 24" w/ Micro Grip <i>Used on 24" Spinergy / Superlight</i>	1	\$272.53
2h	101460	Handrim Natural Fit LT 6 Tab 25" <i>Used on 25" Spinergy/Superlight/Topolino</i>	1	\$234.83
2h	107184	Handrim Natural Fit LT 6 Tab 25" w/ Micro Grip <i>Used on 25" Spinergy / Superlight</i>	1	\$273.67
2i	100950	Handrim Natural Fit LT 6 Tab 26" <i>Used on 26" Spinergy</i>	1	\$234.83
2i	107185	Handrim Natural Fit LT 6 Tab 26" w/ Micro Grip <i>Used on 26" Spinergy</i>	1	\$276.53
3a	200204	Natural Fit Thumb Grip Regular 22"	1	\$144.34
3b	100795	Natural Fit Thumb Grip Regular 24"	1	\$144.34
3c	101462	Natural Fit Thumb Grip Regular 25"	1	\$144.34
3d	200205	Natural Fit Thumb Grip Regular 26"	1	\$144.34
4a	100887	Natural Fit Thumb Grip Super 22"	1	\$189.55
4b	100829	Natural Fit Thumb Grip Super 24"	1	\$189.55
4c	101461	Natural Fit Thumb Grip Super 25"	1	\$189.55
4d	100908	Natural Fit Thumb Grip Super 26"	1	\$189.55
5-7	502169	Handrim Hardware Kit, 6 Rivet Nut, Natural Fit <i>Used with Ki Spoke and MAXX Spoke Wheels</i>	1	\$38.61
8-9	502172	Handrim Hardware Kit, 5 Rivet Nut, Natural Fit <i>Used with 5-Spoke X Core Wheels</i>	1	\$33.10
10	502170	Handrim Hardware Kit, Superlight, EACH <i>Used with Superlight Wheels</i>	1	\$16.55
11-12	502174	Handrim Hardware Kit, 6 Tab, EACH <i>Used with Spinergy and Topolino Wheels</i>	1	\$16.55
5	100698	Handrim Tab	6	\$4.17

Item Number	Part Number	Description	Quantity	Retail
		<i>Used with Ki Spoke and MAXX Spoke wheels. Tabs are included / installed with MAXX Spoke wheels.</i>		
6	100835	10/32 x 1 3/4" Pan Head Screw <i>Used with Ki Spoke and MAXX Spoke Wheels</i>	6	\$1.66
7	100653	Handrim Spacer Spoke <i>Used with Ki Spoke and MAXX Spoke Wheels</i>	6	\$1.10
8	100629	Handrim Spacer Alum Ano Riv <i>Used with 5-Spoke X Core Wheels</i>	5	\$1.16
9	101893	10-32 x 2 PHTH ZN <i>Used with 5-Spoke X Core Wheels</i>	5	\$4.41
10	100536	M5 x 12 BTN SHCS ZC <i>Used with Superlight Wheels</i>	6	\$2.21
11	100669	M5 x 10 BTN SHCS BLKZC w/Patch <i>Used with Spinergy and Topolino Wheels</i>	6	\$1.66
12	100657	M5 Nylock Nut Blk Zn <i>Used with Spinergy and Topolino Wheels</i>	6	\$0.56

Handrim Hardware Chart

Wheel		Handrim Connection				Handrim								
Wheel	Wheel Part No.	Connection Points	Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim		
18" Ki Spoke ¹	200529	3	100698	100653	Aluminum, Plastic Ctd, Projection: 100654	100206		101106						
20" Ki Spoke ¹	200530	6				Natural Fit: 100835	200536		200542					
22" Ki Spoke ¹	200531						100560		100576	100569	200538	200201		
24" Ki Spoke ¹	200532				200349			100577	200547	100793	200202			
25" Ki Spoke ¹	200533				200350		101870	200548	200539	200540				
26" Ki Spoke ¹	200534				200351		100578	200549	100907	101454				
24" Superlight	101159		6			100536	100754	101161	100836		100830	100828		
25" Superlight	101160							101160	101091		101464	101460		
20" 5-Spoke X Core	101961	5		Aluminum, Projection, Natural Fit: 100629 Plastic Ctd: 101756	100724	101897		101963						
22" 5-Spoke X Core	101962					101898		101964						
24" 5-Spoke X Core	100960					100975		100976	200546					
24" 5-Spoke X Core	100960				101893					100768	100769			
20" Spinerger Spox	See Spinerger Spox Page	6			Screw: 100669 Nut: 100657	103125		103179						
22" Spinerger Spox						100827		100808		100889	100888			
24" Spinerger Spox						100766	101161	100615		100830	100828			
25" Spinerger Spox						100767	101160	100765		101464	101460			
26" Spinerger Spox						101477		101148		200200	100950			
22" Spinerger LX	See Spinerger LX Page	6			Screw: 100669 Nut: 100657	100827		100808		100889	100888			
24" Spinerger LX						100766	101161	100615		100830	100828	See Spinerger LX Page		
25" Spinerger LX						100767	101160	100765		101464	101460			
26" Spinerger LX						101477		101148		200200	100950			
18" MAXX Spoke	109285	3	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection: 100654	100206		101106						
20" MAXX Spoke	108244	6				Natural Fit: 100835	200536		200542					
22" MAXX Spoke	105135						100560		100576	100569	200538	200201		
24" MAXX Spoke	105136				200349		101161	100577	200547	100793	200202			
25" MAXX Spoke	107436				200350		101870	200548	200539	200540				
26" MAXX Spoke	107437				200351		100578	200549	100907	101454				

* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 000240 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels or part number 100666 for the screw on 5-Spoke X Core wheels.

¹ Discontinued