


# Clik Backrest Mount - Growth



*Spinergy Back Release Cable based on seat width The Spinergy Back Release Cable length will measure 1" less than the specified seat width. When replacing back release cable, 101991 is needed.*

Item Number	Part Number	Description	Quantity
1a	102556	Spinergy Back Release Cable 8" Wide Black	1
1b	102557	Spinergy Back Release Cable 8" Wide Blue	1
1c	102558	Spinergy Back Release Cable 8" Wide Red	1
1d	102559	Spinergy Back Release Cable 8" Wide Yellow	1

Item Number	Part Number	Description	Quantity
1e	102560	Spinergy Back Release Cable 8" Wide Green	1
1f	102561	Spinergy Back Release Cable 8" Wide Pink	1
1g	102562	Spinergy Back Release Cable 8" Wide Orange	1
1h	102563	Spinergy Back Release Cable 8" Wide White	1
1i	102564	Spinergy Back Release Cable 9" Wide Black	1
1j	102565	Spinergy Back Release Cable 9" Wide Blue	1
1k	102566	Spinergy Back Release Cable 9" Wide Red	1
1l	102567	Spinergy Back Release Cable 9" Wide Yellow	1
1m	102568	Spinergy Back Release Cable 9" Wide Green	1
1n	102569	Spinergy Back Release Cable 9" Wide Pink	1
1o	102570	Spinergy Back Release Cable 9" Wide Orange	1
1p	102571	Spinergy Back Release Cable 9" Wide White	1
1q	102572	Spinergy Back Release Cable 10" Wide Black	1
1r	102573	Spinergy Back Release Cable 10" Wide Blue	1
1s	102574	Spinergy Back Release Cable 10" Wide Red	1
1t	102575	Spinergy Back Release Cable 10" Wide Yellow	1
1u	102576	Spinergy Back Release Cable 10" Wide Green	1
1v	102577	Spinergy Back Release Cable 10" Wide Pink	1
1w	102578	Spinergy Back Release Cable 10" Wide Orange	1
1x	102579	Spinergy Back Release Cable 10" Wide White	1
1y	102580	Spinergy Back Release Cable 11" Wide Black	1
1z	102581	Spinergy Back Release Cable 11" Wide Blue	1
1aa	102582	Spinergy Back Release Cable 11" Wide Red	1
1ab	102583	Spinergy Back Release Cable 11" Wide Yellow	1
1ac	102584	Spinergy Back Release Cable 11" Wide Green	1
1ad	102585	Spinergy Back Release Cable 11" Wide Pink	1
1ae	102586	Spinergy Back Release Cable 11" Wide Orange	1
1af	102587	Spinergy Back Release Cable 11" Wide White	1

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
1ag	200260	Spinergy Back Release Cable 12" Wide Black	1
1ah	200261	Spinergy Back Release Cable 12" Wide Blue	1
1ai	200262	Spinergy Back Release Cable 12" Wide Red	1
1aj	200263	Spinergy Back Release Cable 12" Wide Yellow	1
1ak	200264	Spinergy Back Release Cable 12" Wide Green	1
1al	200265	Spinergy Back Release Cable 12" Wide Pink	1
1am	200266	Spinergy Back Release Cable 12" Wide Orange	1
1an	200267	Spinergy Back Release Cable 12" Wide White	1
1ao	200268	Spinergy Back Release Cable 13" Wide Black	1
1ap	200269	Spinergy Back Release Cable 13" Wide Blue	1
1aq	200270	Spinergy Back Release Cable 13" Wide Red	1
1ar	200271	Spinergy Back Release Cable 13" Wide Yellow	1
1as	200272	Spinergy Back Release Cable 13" Wide Green	1
1at	200273	Spinergy Back Release Cable 13" Wide Pink	1
1au	200274	Spinergy Back Release Cable 13" Wide Orange	1
1av	200275	Spinergy Back Release Cable 13" Wide White	1
1aw	200276	Spinergy Back Release Cable 14" Wide Black	1
1ax	200277	Spinergy Back Release Cable 14" Wide Blue	1
1ay	200278	Spinergy Back Release Cable 14" Wide Red	1
1az	200279	Spinergy Back Release Cable 14" Wide Yellow	1
1ba	200280	Spinergy Back Release Cable 14" Wide Green	1
1bb	200281	Spinergy Back Release Cable 14" Wide Pink	1
1bc	200282	Spinergy Back Release Cable 14" Wide Orange	1
1bd	200283	Spinergy Back Release Cable 14" Wide White	1
1be	200284	Spinergy Back Release Cable 15" Wide Black	1
1bf	200285	Spinergy Back Release Cable 15" Wide Blue	1
1bg	200286	Spinergy Back Release Cable 15" Wide Red	1
1bh	200287	Spinergy Back Release Cable 15" Wide Yellow	1

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
1bi	200288	Spinergy Back Release Cable 15" Wide Green	1
1bj	200289	Spinergy Back Release Cable 15" Wide Pink	1
1bk	200290	Spinergy Back Release Cable 15" Wide Orange	1
1bl	200291	Spinergy Back Release Cable 15" Wide White	1
1bm	200292	Spinergy Back Release Cable 16" Wide Black	1
1bn	200293	Spinergy Back Release Cable 16" Wide Blue	1
1bo	200294	Spinergy Back Release Cable 16" Wide Red	1
1bp	200295	Spinergy Back Release Cable 16" Wide Yellow	1
1bq	200296	Spinergy Back Release Cable 16" Wide Green	1
1br	200297	Spinergy Back Release Cable 16" Wide Pink	1
1bs	200298	Spinergy Back Release Cable 16" Wide Orange	1
1bt	200299	Spinergy Back Release Cable 16" Wide White	1
1bu	200300	Spinergy Back Release Cable 17" Wide Black	1
1bv	200301	Spinergy Back Release Cable 17" Wide Blue	1
1bw	200302	Spinergy Back Release Cable 17" Wide Red	1
1bx	200303	Spinergy Back Release Cable 17" Wide Yellow	1
1by	200304	Spinergy Back Release Cable 17" Wide Green	1
1bz	200305	Spinergy Back Release Cable 17" Wide Pink	1
1ca	200306	Spinergy Back Release Cable 17" Wide Orange	1
1cb	200307	Spinergy Back Release Cable 17" Wide White	1
1cc	200308	Spinergy Back Release Cable 18" Wide Black	1
1cd	200309	Spinergy Back Release Cable 18" Wide Blue	1
1ce	200310	Spinergy Back Release Cable 18" Wide Red	1
1cf	200311	Spinergy Back Release Cable 18" Wide Yellow	1
1cg	200312	Spinergy Back Release Cable 18" Wide Green	1
1ch	200313	Spinergy Back Release Cable 18" Wide Pink	1
1ci	200314	Spinergy Back Release Cable 18" Wide Orange	1
1cj	200315	Spinergy Back Release Cable 18" Wide White	1

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
2	000965	Backrest Cable Release Spring RG	2
3	000909	Backrest Lockout Pin RG	2
4	101991	Spinergy Release Cable/SPOX Adjustment Tool <i>No charge when ordered with Spinergy Release Cable</i>	1