

EHSP-600

AAMI LEVEL-3 – CAT-III TYPE 5B/6B COVERALL



EHSP
Innovation.Speed.Accuracy
for Your Protection



Product Name:	AAMI Level -III CAT-III TYP 5B/6B COVERALL
Product Description:	Breathable, Lightweight and flexible material, particle-tight and spray-tight EN 14126 Coverall.
Manufacturer:	ECEP HAN
Product SKU:	L3/W55GSM -CVR-2020
Brand Name:	EHSP
Product Model:	EHSP 600
Product Color:	White
Private label	Yes Possible, Shall be coordinated with ECEP HAN
Goods in Stock?	Yes
Sizes Available	S, M, L, XL, 2XL, 3XL
Cuff Style	Elasticated
Fastening Type	Zipper Front with Lap Sticker, Elasticated over face, wrists, ankle and waistband.
Height	Varies based on size
Gender	Unisex
Product Material	PP+PE (55 GSM)
Sterility	Non-sterile
Usage	Disposable
User	Adult
Other Features	EN 14126 Standard Defined Fabric protected against infectious diseases, micro - organism such as fluids, aerosols or solid dust particles, this product material is protected, not stipulated by seam technology. EHSP 600 is tight seam sewn, overlapped all seams, elasticated on face, wrist, rear and ankles, the zipper is furnished with metal stopper and flap sticker.

Ecep Han
Trading for good.

ECEP HAN MEDİKAL TEKSTİL SAN VE TİC.LTD.ŞTİ.

Akçaburgaz Mahallesi
Akçaburgaz Caddesi
No:20 AD Suite:28
AKTİM 1 Ticaret ve İş Merkezi
34522-Esenyurt/İSTANBUL/TURKEY



Turkey@ecephan.com
P:+90-212.852 3447
P:+90-212.852 3448
P:+90-212.852 3428

ECEP HAN NORTH AMERICA LLC

8260 E. Gelding Dr. Suite 101
Scottsdale AZ 85260,USA
America@ecephan.com
P:+1-480.656-5519



www.ecephan.com