


(Discontinued 1) Focus CR Stroller Handle Backrest



This Focus CR Stroller Handle Backrest was effective prior to October 2015.

Item Number	Part Number	Description	Quantity
1, 4, 8, 14 - 16	501094	Backrest Hardware Kit Std Focus	1
1, 8, 14 - 16, 18, 19	501406	Backrest Hardware Kit, Std Non Folding, Pair Focus	1
1	001071	Back Plate Std TILT	2
4	001075	Stop Pin Backrest Tilt	1

Item Number	Part Number	Description	Quantity
8	000055	Saddle 1"	4
9	100804	Tube End Plug 1"	1
14	101781	M6 x 50mm Partial Thread BHCS Blk ZC	4
15	100658	M6 Nylock Nut Blk Zn	4
16	100746	M6 Flat Washer Blk Zn	4
18	101815	M8 x 50 BHCS Blk Zn	1
19	100659	M8 Nylock Nut Blk Zn	1