




## (Discontinued 1) Rogue Push Handle



| Size   | Height Range | Length |
|--------|--------------|--------|
| Short  | 11"-14"      | 8.5"   |
| Medium | 13"-17"      | 11.5"  |
| Tall   | 17"-20"      | 15.5"  |



Low Profile  
Push Handle



Standard  
Push Handle

*This Rogue Push Handle is discontinued and was used on serial numbers prior to RG0016612.*

| Item Number      | Part Number | Description                                                                                                                                                           | Quantity |
|------------------|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| 1a, 2, 3a, 4 - 6 | 500779      | Standard Push Handle w/Hardware Assembly<br>Short Rigid<br><i>Use for all fixed heights w/standard push handle and short height adjustable w/standard push handle</i> | 1        |
| 1b, 2, 3a, 4 - 6 | 506730      | Standard Push Handle w/Hardware Assembly<br>Medium Rigid                                                                                                              | 1        |

| Item Number      | Part Number | Description                                                                                                                                                                   | Quantity |
|------------------|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| 1c, 2, 3a, 4 - 6 | 506731      | Standard Push Handle w/Hardware Assembly<br>Tall Rigid                                                                                                                        | 1        |
| 1a, 2, 3b        | 125239      | Low Profile Push Handle Short Rigid<br><i>Use for all fixed heights w/low profile push handle and short height adjustable w/low profile push handle. Replaced part 002758</i> | 1        |
| 1b, 2, 3b        | 506732      | Back Post Assembly Medium with Low Profile<br>Push Handle Rigid                                                                                                               | 1        |
| 1c, 2, 3b        | 506733      | Back Post Assembly Tall with Low Profile Push<br>Handle Rigid                                                                                                                 | 1        |
| 2                | 100638      | Tube End Plug 3/4"                                                                                                                                                            | 1        |
| 3a               | 109432      | Ergo Push Handle Grip 3/4" Standard with<br>Instruction                                                                                                                       | 1        |
| 3b               | 109430      | Ergo Push Handle Grip 3/4" Low Profile with<br>Instruction                                                                                                                    | 1        |
| 4                | 000073      | M5 x 24 Threaded Barrel                                                                                                                                                       | 1        |
| 5                | 100669      | M5 x 10 BTN SHCS BLKZC w/Patch                                                                                                                                                | 2        |
| 6                | 100734      | M5 x 16 PHMS ZC<br><i>Used to mount back upholstery</i>                                                                                                                       | 1        |