

SUPERSTAND

Pediatric Standing System

Multi-Position + Adjustment - Proper Positioning & Comfort



The Art of Standing



ALL IN ONE PRODUCT!

PRIME ENGINEERING

Superstand

MULTI-POSITION STANDER...The Right Approach To Standing

Prone
Upright
Supine

Position Conversions

Conversions from Prone, Upright or Supine are quick and easy. Each Superstand includes a headrest, chin rest and reversible foot system allowing for tool-free position changes.

Fully Adjustable

Support pads, hardware, support straps and chest harness are all adjustable, supplying the best standing position, support and comfort in all configurations.

Tray

Height, depth and angle adjustment allows for optimum positioning of the upper extremities. (Optional Item)

Foot System

Height, width, depth, angle and rotation adjustments allow for best possible foot placement. Entire system is reversible for simple and easy standing position changes.

2 Knee Systems


Adjusts for width, depth and height for optimum positioning and comfort for the lower extremities. (standard)

Optional: Size-Adjustable Kneepads allow you to adjust for various knee sizes and accommodate for braces.

Mobile Base

4 easy rolling casters make the Superstand easy to move and place.

Superstand pictured shows optional tray assembly.



Adjustability and Comfort

Tool-free adjustments allow easy changes to our body support pads, headrest, chinrest, footrest, knee supports, support straps, chest harness and standing angle. Wide range of standing angle adjustment allows for ease of client positioning, adjustment and weight loading. Soft padded positioning supports and cushioned straps, plus complete adjustability, properly position the client to maximize standing time and tolerance.

Growth

Superstand is engineered to grow with your standing needs. **Superstand** fits clients from 26" - 48" and **Superstand Youth** fits clients from 40" - 60" allowing for maximum growth throughout each individual standing program, while maintaining the full positioning capabilities of the **Superstand**. Multiple choices in pad mounting brackets allow for size and growth changes and orthopedic positioning needs.

Features

Superstand Pediatric Multi-Position Stander gives you the right approach to standing - namely *your* approach. **Superstand's** standing range adjustment conveniently operates with a foot pedal allowing for the operators hands to be free. Full adjustability allows you to custom fit **Superstand** to each individual and give them the standing program they need. Growth, Adjustability, Durability, Function, and Ease of Use make the **Superstand** the ideal stander for all your standing needs.

PRIME ENGINEERING 4202 W. SIERRA MADRE AVE., FRESNO, CA 93722 1-800-827-8263 Fax (800) 800-3355
A Division of Axiom Industries Inc. E-mail: Info@primeengineering.com Web Site: www.primeengineering.com

Doc# 70103SS

All specifications subject to change without notice

