

Raz-AP Attendant Propel

MSRP – USA – USD Effective Nov 13, 2023

Dealer / Location:		T:	F:	
·				
PO# / Client:	Contact:	Fmail·		

Raz-AP: Attendant Propel Mobile Shower Commode Chair with height-adjustable stainless-steel frame, Fixed Back, flip-up padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. User weight capacity 407 lbs. For help with justification, please see our <u>Evaluation and Medical Justification Form</u>.

Need help with this order form? Call Raz customer service at (877) 720-5678

1. Standard Raz-AP Frame (18"W) – If you require a special width, depth and / or height see [Sec. 2]

Z10018 Raz-AP w/ 18"W base frame – Standard base frame width, will fit over most toilets	\$1,695
--	---------

2. Custom / Special Frame Widths, Depths and Heights – Adjusta-Back must be selected for all custom widths see [Sec.7]

Z10016	Raz-AP w/ 16"W frame width1 (base chair price – custom upcharge below will be added)	\$1,695					
Z1SW16	Upcharge for custom 16"W frame – will not fit over most toilets	\$495					
Z10020	Raz-AP w/ 20"W frame width ² (base chair price – custom upcharge below will be added)	\$1,695					
Z1SW20	Upcharge for custom 20"W frame	\$495					
Z10022	Raz-AP w/ 22"W frame width (base chair price – custom upcharge below will be added)	\$1,695					
Z1SW22	Upcharge for custom 22"W frame	\$595					
Custom Frame C	onfigurations						
ZLNGF3	Upcharge for custom 3" Longer Frame (3" longer at the front)	\$395					
Z1LOWF4	Upcharge for custom Low Frame (16"-18" STF, 4" lower than standard)	\$395					
Z1APRF3	Upcharge for Recliner Frame (3" longer at the rear for added stability)	\$395					
¹ 18"W Visco-Foan	¹ 18"W Visco-Foam Interface seats can be modified to fit.						
² 18"W seats can b	e selected but must also select Seat Support Kit ZSSKAS [Sec.3]						

DOC#: RSD055r03 Rev. Date: 2023-08-25

3. Seats – Standard – Please note that the Molded and VFI seats, formerly designated 16"D and 19"D, are now designated 17"D and 20"D, which is the center of the IPAS range and the position in which they are shipped from the factory.

Molded Seats (thickness ~2.5"-3")		Visco-Foam Interface Seats (thickness ~3.5")				
ZMS17FRO	18"W x 17"D – Front Access	Standard	ZVFI17BRI	18"W x 17"D – Front Bridged	\$425		
ZMS17BRI	18"W x 17"D – Front Bridged	N/C	ZVFI17FRO	18"W x 17"D – Front Access	\$425		
ZMS20FRO	18"W x 20"D – Front Access	\$135	ZVFI17LTO	18"W x 17"D – Left Access	\$425		
ZMS20BRI	18"W x 20"D – Front Bridged	\$135	ZVFI17RTO	18"W x 17"D – Right Access	\$425		
Specialty Seats	3		ZVFI17RRO	18"W x 17"D – Rear Access	\$425		
ZSOM	Shower-Only Seat – Molded	\$425	ZVFI20BRI	18"W x 20"D – Front Bridged	\$495		
	18"W x 16"-21"D		ZVFI20FRO	18"W x 20"D – Front Access	\$495		
ZEJRM	E&J Replica Seat – Molded – 18"W x 17"D	\$195	ZVFI20LTO	18"W x 20"D – Left Access	\$495		
7510) (-	6005	ZVFI20RTO	18"W x 20"D – Right Access	\$495		
ZEJRV	E&J Replica Seat – Visco Foam – 18"W x 17"D	\$895	ZVFI20RRO	18"W x 20"D – Rear Access	\$495		
Seat Options							
ZSSKAS Seat Support Kit for Raz-AP - These are used to support the sides of seats that are narrower than seat frames on which they are mounted. The reason for a narrower seat on a wider frame is typically to select a Molded Seat and avoid the upcharge of a custom seat. A typical example would be when a 20"W chair is needed but an 18"W standard seat (Molded, VFI, EJR, or Shower Only Seat) is selected.							
ZSS05	Spacers for Front Seat Brackets – 1	L/2" (pr) – for a	ndding dump angl	e (based on 18"D seat)	\$85		
ZSS10	Spacers for Front Seat Brackets – 1	L" (pr) – for add	ding dump angle	(based on 18"D seat)	\$95		

4. Seats - Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. See <u>Custom Seat Assessment Form</u> – Resources tab at razdesigninc.com

ZCSEAT	Custom Dual-Visco Foam Seat (enter specs below) Select if both width and depth are ≤ 21"	\$795
ZCEJR	Custom E&J Replica Seat (enter specs below) Select if both width and depth are ≤ 21"	\$995
ZCSOV	Custom Dual-Visco Foam Shower-Only Seat (no aperture, no drain hole)	\$595
ZCSEATHD	Custom Dual-Visco Foam Seat – HD (enter specs below) Select if width and/or depth is ≥ 22"	\$1,095

Custom Specifications (for above custom seats)

Base Frame Width: __ 16" / 18" / 20" / 22"

Custom Seat Width: __ 16" / 18" / 20" / 22" / 24" (cannot be more than 2" wider than frame width)

Custom Seat Depth: __ 15" / 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" [≥ 21" requires Long Frame (Z1LNGF3)]

Access Opening: __ FR - Front / LT - Left / RT - Right / RR - Rear / BR - Bridged Access Opening Width: __ 2.5" / 3.0" / 3.5" / 4.0" / 4.5" / 5.0" (3.5"W is standard)

Aperture Width: __ 3.0"-8.0" (6.5" recommended) – enter "NA" if aperture is not required

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom Seat width should not be narrower than frame width.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support without ordering wider chair frame.
- For Custom Seats, the Ischial Pelvic Alignment System (IPAS) allows +/- 1" fore / aft adjustment from specified custom seat depth.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accommodation, etc.

5. Fixed Back Frames (18"W only)

This is called a "Fixed" Back because the arm supports are not adjustable in height. Includes molded vinyl back support, fixed-height flip-up padded arm supports, back angle adjustment (90°/97°) with cane-style posts. Back height is measured from top of molded seat to top of back support. Hip Adductors [ZHPG] are not compatible with Fixed Back.

ZFB1818	Fixed Back – 18"W x 18"H (Arm Support Lock Kit (ZALK) [Sec.10] is not included.)	Standard
ZFBU1818	Tension-Adjustable Back Upholstery for Fixed Back	\$185

6. Adjusta-Back Frames - Back frame widths will be the same as the seat frame widths unless otherwise selected in Sec 8

Includes height- and angle-adjustable arms with Armrest Lock Kit (ZALK), tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (90°/97°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. Adjusta-Back – 24"H includes straight posts, on which the upholstery cannot be raised.

ZAB 15*	Adjusta-Back ¹ – 15"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	\$225
ZAB 17*	Adjusta-Back ¹ – 17"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	\$225
ZAB 20*	Adjusta-Back ¹ – 20"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	\$295
ZAB 24*	Adjusta-Back ¹ – 24"H	* insert back frame width to complete the code (16" / 18" / 20" / 22" / 24")	\$345

Reclining Backs – 30° (must select Reclining Frame option Z1APRF3 in Sec 2) Reclining back frame width must be the same as the seat frame width. 2-Point Arms (Z2PLAS / Z2PSAS) are included with Reclining Backs.

ZRBA 20*	Reclining Adjusta-Back with 20"H back (anti-tippers included) * insert back frame width to complete the code (16" / 18" / 20" / 22")	\$1,095
ZRBA 24*	Reclining Adjusta-Back with 24"H back (anti-tippers included) * insert back frame width to complete the code (16" / 18" / 20" / 22")	\$1,195
ZDTL	Dual Triggers – both on left side	N/C
ZDTR	Dual Triggers – both on right side	N/C

 $^{^{1}}$ Dual Tilt Triggers option positions both triggers on one side, allowing for single-handed activation.

7. Custom Back Frame Width Narrower Than Chair Frame Width – Must choose back frame [Sec.7] (configurations below are not available for recliner backs)

NOTE: For back	NOTE: For back frame width variances on recliner frames, please request custom quotation.								
ZNB	_	ustom – Back Frame 2"	\$345	Chair Frame Width	16"	18"	20"	22"	
	N	larrower than chair frame		Back Frame Width	14" ¹	16" ²	18" ³	20'' ⁴	
* Enter chair and back frame widths to complete code.							d b		
eg. ZNB2018 sp	pecifies a 2	0"W chair frame with							
an 18"W back	frame								

¹ Requires Custom Seat

² Compatible with all 18"W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.

³ Custom Seat recommended. Also compatible with all 18"W seats with Seat Support Kit Z4SSKAS.

⁴ Custom Seat recommended.

8. Back Support Options – These back supports replace the molded vinyl upholstery on a Fixed Back or tension-adjustable upholstery on an Adjusta-Back. Symphony and Harmony Backs are recommended for infection-control environments.

ZSB14	14"W Symphony Back ¹ – recessed solid back support	\$295	
ZSB16	16"W Symphony Back ¹ – recessed solid back support	\$295	
ZSB18	18"W Symphony Back ¹ – recessed solid back support	\$295	
ZSB20	20"W Symphony Back ¹ – recessed solid back support	\$395	
ZSB22	22"W Symphony Back ¹ – recessed solid back support	\$465	
ZHB	Harmony Back $^{\!\!\!1,2}$ (contoured solid back – 18"W x 13"H) Fits 18"W back frames only Fits 17" / 20" back heights only	\$295	
¹ Does not fit 15"H backs			
² Not recomm	nended for 24"H backs or reclining backs.		

9. Arm Supports and Accessories

	ZARML	Long Arm Supports (pr) – 13"L pad	Std		ZPAML	3	\$165
	ZARMS	Short Arm Supports ¹ (pr) – 10"L pad	N/C			(for ZFPLL/ZTRL)	
	ZALK	Arm Support Locks (pr) for Fixed Back (ZFB1818)	\$75		ZPAMR	Pivoting Arm Mount – Right (for ZFPLR / ZTRR)	\$165
	Z2PLAS	2-Point Arm Supports ² (pr) with 13"L pads (std with reclining backs)	\$345		ZASK05	Arm Support Spacer Kit – 1/2" per side for all Adjusta-Backs (pr)	\$60
	Z2PSAS	2-Point Arm Supports ² (pr) with 10"L pads (std with reclining backs)	\$345		ZASK10	Arm Support Spacer Kit – 1" per side for all Adjusta-Backs (pr)	\$115
	ZFAPLL		\$90		ZPHGL	Pivoting Hand Grips (pr) for 13"L arms ⁴	\$395
	ZFAPLL	Long Molded Flat Armpad ³ – Left (14"L x 3.25"W – fits ZARM13)	\$90		ZPHGS	Pivoting Hand Grips (pr) for 10"L arms ⁴	\$395
	ZFAPLR	Long Molded Flat Armpad ³ – Right (14"L x 3.25"W – fits ZARM13)	\$90		ZAPSB*	* Enter frame width to complete code	\$225
	ZFAPSL	Short Molded Flat Armpad ³ – Left $(10^{\circ}L \times 3^{\circ}W - \text{fits ZARM10})$	\$90		ZAPSB22	/co	\$245
	ZFAPSR	Short Molded Flat Armpad ³ – Right $(10^{\circ}L \times 3^{\circ}W - \text{fits ZARM10})$	\$90			for a 22"W Back Frame	
	ZTRL	Molded Arm Trough ³ – Left (18.5″L x 4.5″W)	\$110				
	ZTRR	Molded Arm Trough ³ – Right (18.5"L x 4.5"W)	\$110				
¹ Not compatible with ZFAPL / ZFAPS / ZTR / ZPAM							
	² For chairs with same back width and frame width only					combined with 2 Point Arms requires a custom	
3	May require Z	PAM / ZASK05 / ZASK10 to flip up, depending on position of the Flat Arm Pad or Arm Trough		Ι΄	ocharge. Not compatible	e with ZFAPL / ZFAPS / ZTR	

10. Foot and Leg Supports

_						
	ZFS	Foot Supports 16"-20.5" (pr)	Std	ZELRL	Elevating Leg Support ³ – Left	\$495
	ZFS16	Foot Supports for 16"W frames (pr)	N/C		(includes ZADML)	
		(16"-20.5")		ZELRR	Elevating Leg Support ³ – Right	\$495
	ZFET2	Footrest Extension Tubes ¹ (pr) – 2" longer	\$50		(includes ZADMR)	
	ZMFX8	MFX8 Foot Supports 8"-12" (pr)	\$40	ZFOR	Forward Offset Foot Support Receiver –	\$95
	ZMFX12	MFX12 Foot Supports 12"-16.5" (pr)	\$70		(ea) Moves foot support forward 2.5" to	
	ZADM	Adj Angle/Depth Molded Footplates ² (pr)	\$120		accomodate deeper seats and reduced knee flexion range. Specify 1 or 2	
	ZFET2A	Footrest Extension Tubes (pr)	\$50	ZRLSL	Universal Residual Limb Support – Left	\$195
		for ZADM – 2" longer		ZRLSR	Universal Residual Limb Support – Right	\$195
	ZADABS16	Angle/Depth Adj ABS Footplates For 16"W frames (pr)	\$195			
	ZVS*	V-Style Foot Support (1-piece) 16"-20.5" * Enter frame width to complete code (16" / 18")	\$295			
	ZVS*	V-Style Foot Support (1-piece) 16"-20.5" * Enter frame width to complete code (20" / 22")	\$345			
	ZFBAS	Flip-Back Foot Support (1-piece) 16"-20.5"	\$295			
	Not compatible	e with ZMFX12 / ZADM / ZVS / ZFBT			patible with frame widths \geq 18" and adjusted sea or \leq 21" w/3" longer seat frames (ZLNGF3)	at

11. Foot and Leg Support Accessories

ZFSLK	Foot Support Lock (pr)	\$95	ZNC	Neoprene Covers for Standard	\$90
ZCS *	Calf Strap, Snap-On – Bodypoint	\$105		Footplates (pr)	
	Aeromesh * Enter frame width to complete code (16" / 18" /20" / 22")		ZNCA	Neoprene Covers for Angle/Depth Adjustable Footplates	\$90
ZCSIC*	* Enter frame width to complete code	\$95	ZNC16	Neoprene Covers for Narrow Molded Footplates – 16"W frame (pr)	\$90
	(16" / 18" /20")		ZNCABS16	Neoprene Cover for Narrow Molded ABS	\$90
ZCSIC22	Calf Strap, Snap-On – Infection Control	\$125		Angle/Depth Adjustable Footplates	
	For a 22"W frame		ZNCVS	Neoprene Cover for V-Style Foot Support	\$90
ZCPL	Calf Panel – Bodypoint Aeromesh ¹ Fits 16" / 18" /20"W frames	\$225	ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$90
ZHS1820	H-Strap – for 18" / 20"W frames – hook and loop attachment	\$90	ZLMOR	Lateral / Medial Offset Foot Support Receiver (ea) – moves foot support 2"	
ZHS2224	H-Strap – for 22"W frames	\$110		laterally or medially. Specify 1 or 2	
	 hook and loop attachment 		ZFOR	Forward Offset Foot Support Receiver	\$95
ZAHL	Adjustable Heel Loops (pr) (not available for ZVS / ZFBAS)	\$60		(ea) – Moves foot support forward 2.5" to accomodate compromised knee flexion.	·
¹ Not compatible with ZMFX12 / ZFBAS / ZVS			Specify 1 or 2		

12. Positioning Pads – Pad Sizes: Small 5.5" x 3.5" / Large – 6.75" x 4.3".

ZHSAS2	Head Support with small pad	\$495			
ZHSAS2		\$495			
	Head Support with large pad				
ZOHIP	Head Support with large pad				
	Lateral Offset Head Support Interface Plate – up t	o 22° of adjustment \$75			
	cic Supports – Pads on Lateral Thoracic Supset, specify ZLAEP extension plate.	oports can be adjusted medially up to 3" per side. For up to			
ZLTS1L	Lateral Thoracic Support w/ Small Pad – Left	\$225			
ZLTS1R	Lateral Thoracic Support w/ Small Pad – Right	\$225			
ZLTS2L	Lateral Thoracic Support w/ Large Pad – Left	\$225			
ZLTS2R	Lateral Thoracic Support w/ Large Pad – Right	\$225			
Leg Adductor	s				
ZLGAS1L	Leg Adductor w/ Small Pad – Left	\$225			
ZLGAS1R	Leg Adductor w/ Small Pad – Right	\$225			
ZLGAS2L	Leg Adductor w/ Large Pad – Left	\$225			
ZLGAS2R	Leg Adductor w/ Large Pad – Right	\$225			
Hip Adductors	S				
	Hip Adductor w/ Small Pad – Left	\$225			
ZHP1R	Hip Adductor w/ Small Pad – Right \$225				
ZHP2L	Hip Adductor w/ Large Pad – Left \$225				
ZHP2R	Hip Adductor w/ Large Pad – Right \$225				
Hip Adductors F	or Reclining Backs				
ZHPR1L	Hip Adductor for Reclining Back w/ Small Pad – Left \$225				
ZHPR1R	Hip Adductor for Reclining Back w/ Small Pad – Right \$225				
ZHPR2L	Hip Adductor for Reclining Back w/ Large Pad – Left \$225				
ZHPR2R	Hip Adductor for Reclining Back w/ Large Pad – Right \$225				
	s for Lateral Thoracic Supports and Hip / Leg Ado ow for an additional 1.5" medially and can be use	luctors d on Lateral Thoracic Support, Hip and Leg Adductors			
ZLAEP	Lateral/Adductor Extension Plate (ea)	\$28			
Lateral Thoracic	Support Left Right				
Leg Adductor	Left Right				
Hip Adductor	Left Right				
Total Quantity _					
Medial Thigh Su	ipports (Pommel)				
ZMTSAS1	Medial Thigh Support w/ Small Pads \$595				
ZMTSAS2	ZMTSAS2 Medial Thigh Support w/ Large Pads \$595				

13. Pelvic and Chest Belts - A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts		_				Che	est Belts					
ZPBS	Pelvic Belt - 2-piece w/	•	rame widths 14- ase buckle)	-16",	\$50		ZCBS		Belt ² – Small (bac e w/ hook and lo		l"-16",	\$45
ZPBM			n (frame widths -release buckle)		\$50		ZCBM		Belt ² – Medium (piece w/ hook ar			\$45
ZPBL	Pelvic Belt - 2-piece w/	• .	rame widths 18' ase buckle)	"-22",	\$50		ZCBL		Belt ² – Large (bade w/ hook and lo		8"-22",	\$45
ZPBXL		•	seat frame widt side-release buc		\$65		ZCBXL	Chest	Belt – XLarge ² (ba 24", 1-piece w/ ho	ack frame w		\$55
ZPBICS			on Control, Smal	II	\$75			closur	•	50K ana 100	P	
	(frame widew/ side-rele		•				ZCBICS		Belt ^{1,2} – Infection e w/ side-release	,		\$75
ZPBICM	retrie Beit intection control, Wedian				widths 14-16")							
	(frame widew/ side-rele		•				ZCBICM	Mediu	Chest Belt ^{1,2} – Infection Control, Medium 1-piece w/ side-release buckle		uckle	\$75
ZPBICL	Pelvic Belt ¹ – Infection Control, Large (frame widths 18"-22", 2-piece w/ side-release buckle)		е	\$85			•	widths 16"-18")				
					ZCBICL	1-piece w/ side-release buckle (back			\$85			
ZPBICXL			on Control, XLar	ge	\$95		ZCDICVI		s 18"-22")			ĆOF
	w/ side-rele		2"-24", 2-piece le)				ZCBICXL Chest Belt ^{1,2} – Infecti 1-piece w/ side-relea		e w/ side-release	buckle (ba		\$95
ZPBBPS	Pelvic Belt -	– Bodypo	int Aeromesh, S	mall	\$245		760006		widths 22"/ 24")		C II	ć275
ZPBBPM	Pelvic Belt -	– Bodypo	int Aeromesh, N	∕led	\$215		ZCBBPS		Belt – Bodypoint	·		\$275
ZPBBPL	Pelvic Belt -	– Bodypo	int Aeromesh, L	arge	\$215		ZCBBPM		Belt – Bodypoint			\$245
ZPBBPXL	Pelvic Belt -	– Bodypo	int Aeromesh, X	Large	\$175		ZCBBPL		Belt – Bodypoint		_	\$245
						<u> </u>	ZCBBPXL		Belt – Bodypoint	1		\$195
Recommended sizes for	Pelvic belt	ZPBBPS	N/A	ZPBBPN		'W se frai	me	ZPBBPL	18"/20"W base frame	ZPBBPXL	22"/24"V base fran	
Bodypoint belts	Chest belt	ZCBBPS	14" back frame	ZCBBPN		'W ck fran	me	ZCBBPL	18"/20"W back frame	ZCBBPXL	22"/24"V back fram	
1												

 $^{^{1}}$ Non-absorbent belts – recommended for infection-control environments

14. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$67			
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$125			
ZSC17	Splash Curtain for 17"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat $1/2$ ").	\$325			
ZSC20	Splash Curtain for 20"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat $1/2$ ").	\$345			
ZSCFR	Splash Curtain for 17" / 20"D (at front of aperture only, not for Front Open Seats)	\$345			
ZSCC	Custom Splash Curtain for Custom Seats – Raises seat 1/2"	\$345			
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$345			
ZBUDBH	Bard urine drainage bag holder (available for all frame widths)	\$155			
¹ Bridged Front Se	Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for pricing.				

² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.

15. Other

ZDCLA			ZPSK1	Commode Pan Spacer Kit ² – 1" lower	\$200		
caster (facilitates steering)		\$245	ZPSK2	Commode Pan Spacer Kit ² – 2" lower	\$215		
ZATAS	7 that Tippers (pr/7 that Tippers are sized		ZPSK3	Commode Pan Spacer Kit ² – 3" lower	\$230		
	according to caster diamter ordered. Rear Front		ZPSK4	Commode Pan Spacer Kit ² – 4" lower	\$245		
ZCUAP4	4" Caster Upgrade (4x)	\$120					
ZCUAP6	6" Caster Upgrade (4x)	\$350	² Commode Pan Spacer Kits provide extra clearance for self-care with				
¹ Rear Anti-Tip	¹ Rear Anti-Tippers are included with Reclining Back Option [Sec.7]			Increases height required to clear toilet.	are with		

16. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z907	Omit Foot Supports - if Foot / Leg	- \$45
				Support is not required	

Some assembly of Raz products is required.

17. Notes:		

407 lbs

Specifications

Weight Capacity

Raz-AP (Attendant Propel)

Weight Supusity	107 150
Overall Width	Seat width plus 4.5" (with Standard Arm Support Pads)
Overall Length (foot supports in middle position)	36.5" (with rear caster trailing), 35" (to rear canes)
Seat Depth	16" or 21" depending on seat/cushion length
Toilet Clearance Height (under seat) - with pan brackets - without pan brackets	15.25"-19.25" (1" increments) 16"-20" (1" increments)
Toilet Clearance Width	16" for standard 18"W seat frame. For other widths, seat frame width minus 2"
Arm Support Height - above molded seat (1/2" inc	rements)
 18"H Fixed Back - standard 15"H Adjusta-Back 17"H Adjusta-Back 20"H Adjusta-Back 24"H Adjusta-Back 20"H Reclining Adjusta-Back 24"H Reclining Adjusta-Back 	8.5" 5.5"-8.5" 6.5"-9.5" (1/2" increments) 7.5"-10.5" (1/2" increments) 7.5"-10.5" (1/2" increments) 7.5"-10.5" (1/2" increments) 7.5"-10.5" (1/2" increments)
Seat-to-Floor Height (with molded seat)	20"-24"
Product Weight (with 18"W seat)	38 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)