



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

# R82 Crocodile

Company:	Ship To:		
Account Number:	Date:	Address:	
Quote:	PO Number:	City:	State:
ATP:	Zip:	Phone:	
Marked For:	Email/Fax:		
Associated Quote/Order#:		Effective: 1/20/2025	

## Crocodile

	0	1	2	3
Width between handles	4 - 9¾"	7½ - 13¾"	9½ - 16½"	13 - 19"
Height from floor to handles	11¾ - 24"	16½ - 28¾" 15 - 21¼"*	21¾ - 33"	29½ - 41¼"
Max. user height	41¼"	45¼"	59"	70¾"
Max. user mass	66 lbs	66 lbs	99 lbs	176 lbs
Length	25¼"	28¼"	33"	34¼"
Overall width	22¾"	24¾"	27¼"	27¼"
Folded length	25½"	28½"	33¾"	39¼"
Folded width	22¾"	24½"	27¼"	27¼"
Folded height	11¾"	12¾"	10½"	11½"
Minimum turning radius	19"	19¾"	22"	45"
Mass device	13.4 lbs	15.7 lbs	15.7 lbs	26.7 lbs

\*Mounted with the Starter Kit accessories. Further measurements at our website.

- Width and height adjustable handles
- Easy to fold
- Lightweight durable aluminium frame – Sizes 0, 1 and 2
- Stable steel and aluminium frame – Size 3
- Upholstery and cover for top frame
- Solid front and rear wheels
- Swivel front wheels with directional stabilizer
- Tension adjustable rear wheels
- On/off anti-reverse wheel stop
- Padded Rear Stabilizer Bar w/ Tool Kit



HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.  
Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.  
**DCN1197.2 Pricing and specifications subject to change without notice.**

## Crocodile Frame

### Crocodile Base

#### Frame

Crocodile, White, Size 0 (86800)	\$733
Crocodile, Yellow, Size 1 (86801)	\$808
Crocodile, Orange, Size 2 (86802)	\$808
Crocodile, Metallic Black, Size 2 (86852)	\$808
Crocodile, Anthracite Grey, Size 3 (86803)	\$1333

### Mounting Kits



Crocodile Starter Kit, for Crocodile size 1 (86801-0)	\$513
<i>Used for Size 1</i>	

### Anti-Tippers (Pair)



#### Anti-Tip Device, Swing-away

Anti-tip Swing-away, Pair, for Size 0-2 (86815)	\$175
Anti-tip Swing-away, Pair, for Size 3 (86815-3S)	\$175



#### Anti-Tip Device, Fixed

Anti-tip Fixed, Pair, for Size 1-2 (86815-F)	\$84
Anti-tip Fixed, Pair, for Size 3 (86815-3F)	\$84

## Arm and Forearm Options

### Arm Supports



#### Forearm Support with Handgrip

Forearm Support Pair - for Sizes 0-1 (86858-1)	\$472
Forearm Support Pair - for Size 2 (86859-1)	\$472
Forearm Support Pair - for Size 3 (86860-1)	\$472



#### Padded Forearm Positioner (Each)

*For use with Forearm Support with Handgrip. For qty 2 selection both options below.*

Padded Forearm Positioner Strap (869370)	\$75
Padded Forearm Positioner Strap (869370)	\$75



#### Handle Knobs

Handle Knobs, Pair (86861)	\$56
----------------------------	------

## Handles and Grip Bars

### Grip Bar



#### Grip Bar For Crocodile

Grip Bar, for Size 0 (86879-0)	\$141
Grip Bar, for Size 1 (86879-1)	\$163
Grip Bar, for Size 2 (86879-2)	\$196

### Handles and Push Braces



#### Attendant Guide Bar

Attendant Handle (86866-1)	\$249
----------------------------	-------

#### Angle Adjustable Grip Bar

Angle Adjustable Grip Bar, Size 0 (86878-0)	\$235
Angle Adjustable Grip Bar, Size 1 (86878-1)	\$275
Angle Adjustable Grip Bar, Size 2 (86878-2)	\$305

## Brakes


### Brake

*Hand brake eliminates the function of anti-reverse wheel stop. If ordering hand brake with anti-reverse wheel stop function, a custom order required.*



#### Hand Brake

Hand Brake Kit, for Size 1-2 (86834)	\$262
Hand Brake Kit, for Size 3 (86873)	\$262



#### Deceleration Hand Brake Kit

Deceleration Hand Brake Kit, Sizes 1-2 (86834-DEC)	\$262
Deceleration Hand Brake Kit, Size 3 (86873-DEC)	\$262

### Holder for Hand Brake

*Required with use of Hand brake and Forearm support. Allows for mounting of handbrake to Forearm support.*

#### Holder for Hand Brake

Small, for Size 1 (86874-1)	\$367
Medium, for Size 2 (86874-2)	\$414
Large, for Size 3 (86874-3)	\$445

## Back and Pelvic and Hip Supports

### Back and Pelvic Supports



#### Surrounding Supports

Surrounding Support for Size 0 (868610)	\$535
---	-------



#### Back Support

Back support - Adjustable Padded for Size 3 (86886)	\$552
---	-------



#### Pelvic Support With Laterals

Pelvic Support With Laterals, for Sizes 0-2 (86808)	\$287
Pelvic Support With Laterals, For Sizes 1-2 (86809)	\$287

#### Abdominal Support Belt

*For size 3 back support only*

Abdominal Support Belt, for Size 3 (85490-2)	\$134
--	-------



#### Lateral Strap

Lateral Strap (905011)	\$56
------------------------	------



#### Mounting Hardware

*Required for use of pelvic and back support on size 3*

Cross Bar for Mounting Back/Pelvic Supports, for Size 3 (86880)	\$107
---	-------

### Hip Support Pair



#### Hip Support Pair

Hip Supports, Pair (86819)	\$257
----------------------------	-------



Phone (877) 227-0097 | Fax (704) 882-0751  
 sales.us@etac.com

# R82 Crocodile

## Seating Options

### Solid Flip Down Seat



#### Solid Flip Down Seat

White, Size 0 (86820-01)	\$186
Yellow, Size 1 (86820)	\$186
<i>Not compatible with Crocodile starter kit</i>	
Orange, Size 2 (86821)	\$186
Metallic Black, Size 2 (86823)	\$186
Anthracite Grey, Size 3 (86863)	\$186

### Seat Slings



Seat Sling, for Size 1 (86824)	\$147
Seat Sling, for Size 2 (86825)	\$147



*Requires pelvic support with lateral and lateral strap or surrounding support*

Lightweight Sling Seat, for Size 0-1 (86824)	\$147
Lightweight Sling Seat, for Size 2 (86825-1)	\$147

## Accessories

### Bags and Baskets



#### Bags and Baskets

Storage Bag for Sizes 0-2 (86839)	\$78
Storage Bag, for Size 3 (86843)	\$230