



# Raz-AP

## Attendant Propel Rehab Shower Commode

The Raz-AP (Attendant-Propel) rehab shower commode chair comes standard with 4" (10 cm) of height adjustment, tension-adjustable back upholstery.



FEATURES

Extensive range of options available

Flip-up padded arms

Exclusive to Raz, Ischial Pelvic Alignment System (IPAS), provides 2" (5 cm) of fore/aft or 1" (2.5 cm) of lateral seat adjustment

Stainless steel frame with 4" (10 cm) of height adjustment



Tension-adjustable back upholstery

Adjustable back angle

Dual-locking (swivel and roll) casters

**ALSO AVAILABLE**

- Long frame – 3" (8 cm) longer for tall clients
- Low frame - 4" lower seat-to-floor height
- 16" and 20" (40 and 51 cm) widths
- Recliner frame

| SPECIFICATIONS                                   | Raz-AP (Attendant Propelled)   |
|--|--|
| WEIGHT CAPACITY                                  | 400 lbs (181 kg)   |
| FRAME WIDTH (seat and back frame)                | 18" (46 cm)  |
| OVERALL WIDTH (with standard armpad)             | 22.5" (57 cm)  |
| OVERALL LENGTH (with footrest)                   | 33.75" (86 cm)   |
| SEAT WIDTH (commode seat)                        | 18" (46 cm)  |
| SEAT DEPTH (front of back post to front of seat) | 16" or 19" (40 or 48 cm)   |
| CLEARANCE HEIGHT UNDER SEAT                      | 15.5" - 19.5" (39 – 50 cm) with pan brackets<br>16.25" - 20.25" (41 - 51 cm) without pan bracket |
| ARMREST HEIGHT ABOVE SEAT (17" Adjusta-Back)     | 6.5" - 9.5" in ½" increments (17 - 24 cm)  |
| SEAT-TO-FLOOR HEIGHT                             | 20" - 24" in 1" increments (51 - 61 cm)  |
| PRODUCT WEIGHT (with footrests)                  | 38 lbs (17.25 kg)  |
| FRAME AND HARDWARE MATERIAL                      | Stainless Steel (Lifetime warranty against corrosion)  |

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.  
Specifications, and options are subject to change without prior notice.

MANUFACTURED BY:

RAZ Design Inc. 135 Railside Road, Toronto, Ontario Canada, M3A 1B7 T: (877) 720-5678 E: info@razdesigninc.com