


Axiom DEX Hardware



New DEX Hardware (Effective 5/30/2023) Removable Adjustable Hardware (300lb Max User Weight)

Item Number	Part Number	Description	Quantity	Retail
1-14	119829	DEX Hardware Retrofit Kit, Axiom Back	1	\$340.00
1, 2, 10, 11a	119831	DEX Hardware Housing and Hardware Assembly, Left, Axiom Back	1	\$61.50
1, 2, 10, 11b	119832	DEX Hardware Housing and Hardware Assembly, Right, Axiom Back	1	\$61.50
1, 2, 10, 11a, 11b	119830	DEX Hardware Housing and Hardware Assembly, PAIR, Axiom Back	1	\$123.00
7-9, 12, 13, 14a, 14b	119833	DEX Clamp and Post Assembly, Hardware Included, Each, Axiom Back	1	\$78.50
7-9, 12, 13, 14a, 14b	119834	DEX Clamp and Post Assembly, Hardware Included, PAIR, Axiom Back	1	\$156.50
1	101823	M6 x 20 HHCS Blk Zn	2	\$1.13
2	101853	M6 Washer Blk Zn	2	\$0.58

Item Number	Part Number	Description	Quantity	Retail
3	100746	M6 Flat Washer Blk Zn	2	\$1.71
4	100558	M6 NYL INS Domed Cap L/N ZC	2	\$1.13
7	005476	Nut, Dual, SR.25, M6, .60SP	1	\$38.61
8	005472	Axiom Back Mount Clamp	2	\$33.09
9	005469	Axiom Back Mount Post	1	\$33.09
12	101835	M6 Lock Washer Blk Zn	2	\$0.58
13	003810	Washer, Spher, 10mm	2	\$11.03
14a	100860	M6 x 35 HHCS P/THD ZC <i>Used with 1" and 1-1/8" Back Canes</i>	2	\$1.71
14b	101824	M6 x 30 HHCS Blk Zn <i>Used with 3/4" and 7/8" Back Canes</i>	2	\$2.28