



Breezy 600

EIB60, January 2022

ACCOUNT INFORMATION

Account No: _____ Drop Ship/Ship To: _____
 Date: _____ Name: _____
 P.O. #: _____ Address: _____
 Buyer: _____ City: _____ State: _____
 Marked For: _____ Zip Code: _____ Phone Number: _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

1. FRAME

BREEZY 600

Hangers not included in base price

- 086WM20 K0004 Breezy 600 \$1,980.00

2. SIZE / SEATING

Seat to Floor Height Range

- Standard (15" - 20 1/2", see matrix) Std
 Ultra Hemi (13" - 18", see matrix) NC

Seat Width

- 086W5 14" Std
 086W7 16" Std
 086W9 18" Std
 086W42 E2201 20" \$151.00

Seat Depth

- 086FD7 16" Std
 086FD9 18" NC

Seat Height

- 086HT1 13" \$102.00
 086HT2 14" \$102.00
 086HT3 14 1/2" \$102.00
 086HT4 15" \$102.00
 086HT5 15 1/2" \$102.00
 086HT6 16" \$102.00
 086HT7 16 1/2" \$102.00
 086HT8 17" NC
 086HT9 17.5" NC
 086HT10 18" NC
 086HT11 18.5" NC
 086HT12 19" ¹² NC
 086HT13 19.5" ¹² NC
 086HT14 20" ¹² NC
 086HT16 20.5" ¹² NC

Seat Sling

- 086SL3 Padded Seat Sling \$15.00
 086SL4 Non-padded Seat Sling ² Std
 086SL6 Omit Seat Sling -\$15.00

Cushion

To order any JAY cushions on-chair, please see JAY price list.

- 086JC01 JAY Basic Price List
 086JC02 JAY Easy Price List
 086JC03 JAY Soft Combi P Price List
 086JC04 JAY J2 Price List
 086JC05 JAY J2 Deep Contour Price List
 086JC06 JAY Lite Price List
 086JC07 JAY J3 Price List
 086CU1 K0669 2" Cushion \$161.00

3. CASTERS

CASTERS

- 086CW5 5" Polyurethane NC
 086CW32 8" Polyurethane ¹⁶ Std
 086CW21 8" Pneumatic ⁶ \$92.00

Caster Options

- 086CS9 1-1/2" Fork Stem (Double Bearing) \$38.00
 086CP13 Double Bearing Caster Plate \$38.00

Footnotes:

¹ Not available when using Standard frame.

² For use with cushion.

- ³ Only available with standard hemi extension tube.
⁴ Not available with Ultra Hemi frame.
⁵ Not available with 8" pneumatic casters.
⁶ Not available with 20" rear wheels when using Standard frame.
⁷ Only available with 8" casters.
⁸ Not available with 14" seat width.
⁹ Pull-to-lock wheel locks required.
¹⁰ Not available with composite footplates.
¹¹ Not available with articulating legrest.
¹² Must use 1-1/2" stem.
¹³ Not available with ELR.
¹⁴ Only available with ELR
¹⁵ To order a J3 Back please select Omit Back Upholstery.
Please note: The J3 Back requires back canes or push handles.
¹⁶ Not available with 70° hangers due to interference

4. REAR WHEELS

Rear Wheel Size		
<input type="checkbox"/>	086WS4	20" NC
<input type="checkbox"/>	086WS5	22" NC
<input type="checkbox"/>	086WS6	24" Std

Rear Wheel Tires		
<input type="checkbox"/>	086RT4	Low-profile Polyurethane ^{17,21} Std
<input type="checkbox"/>	086RT29	E2218 Full-profile Polyurethane ^{18,21} \$54.00
<input type="checkbox"/>	086RT35	E2211 E2212 Pneumatic ²¹ \$177.00
<input type="checkbox"/>	086RT33	E2211 E2213 Pneumatic w/Airless Insert ^{21,22} \$230.00

Hand Rims		
<input type="checkbox"/>	086HR5	Composite ²¹ Std
<input type="checkbox"/>	086HR2	Aluminum - painted ²¹ NC
<input type="checkbox"/>	086HR1	Aluminum - anodized ²¹ \$0.00
<input type="checkbox"/>	086HR4	Plastic Coated ²¹ \$118.00
<input type="checkbox"/>	086HR7	Omit Hand Rims ²¹ NC

Axle		
<input type="checkbox"/>	086AL4	Threaded Offset Axles Std
<input type="checkbox"/>	086AX43	Quick-Release Rear Axles \$102.00
<input type="checkbox"/>	086OA2	E0958 One-Arm Drive ^{23,24} \$913.00
<input type="checkbox"/>	086HM2	Left Mount
<input type="checkbox"/>	086HM1	Right Mount

Anti-Tip Tubes		
<input type="checkbox"/>	086AT1	E0971 Rear Anti-Tip Tubes \$161.00

Projections ^{19, 21}		
<input type="checkbox"/>	086HR8	Oblique - 8 \$349.00

Wheel Locks		
<input type="checkbox"/>	086WL1	High-Push Breezy-Style Std
<input type="checkbox"/>	086WL2	High-Pull Breezy-Style NC
<input type="checkbox"/>	086WL22	High-Push Quickie-Style \$43.00
<input type="checkbox"/>	086WL21	High-Pull Quickie-Style \$43.00

Wheel Lock Options		
<input type="checkbox"/>	086EH1	E0961 6" Ext. Handles with Bungee \$81.00
<input type="checkbox"/>	086GA1	E0974 Grade Aids ^{20,25} \$161.00

See Sunrise Seating Solutions on reverse side.

TOTAL CHAIR PRICE: \$ _____

FOOTNOTES:

- ¹⁷ Not available with Projection Hand Rims.
¹⁸ "No Charge" option when used with Projection Hand Rims.
¹⁹ Not available with Low-profile polyurethane tires. Full-profile Polyurethane tires available as "No Charge" option. Requires 24" rear wheels.
²⁰ Not available with Polyurethane tires.
²¹ Not available with 24" Mag SLT
²² Must order Pneumatic tire with Flat Free Insert
²³ Not available with Ultra Hemi frame.
²⁴ Not available with Quick-Release Rear Axles.
²⁵ Grade Aids require Quickie style wheel locks

5. FOOTREST

Hanger Type		
<input type="checkbox"/>	086H4	60° Swing-Away Hemi ^{1,3} \$151.00
<input type="checkbox"/>	086H131	60° Swing-Away ⁴ \$151.00
<input type="checkbox"/>	086H2	70° Swing-Away ^{4,5} \$151.00
<input type="checkbox"/>	086H24	K0053 Articulating Legrests ^{4,6,7,8,9} \$397.00
<input type="checkbox"/>	086H14	E0990 Elevating Legrests ^{4,10} \$338.00

Footrest Extension Tubes		
<input type="checkbox"/>	086E1	Short (13"-15 1/2") ^{4,11} NC
<input type="checkbox"/>	086E2	Medium (15 1/2" - 18") ^{4,6} Std
<input type="checkbox"/>	086E44	Long (18" - 20 1/2") ^{4,6,12} \$38.00
<input type="checkbox"/>	086E21	Omit Ext. Tubes & Footplates NC

Footplates		
<input type="checkbox"/>	086F1	E0951 Composite w/Heel Loops ¹³ Std
<input type="checkbox"/>	086F5	Aluminum ¹⁴ NC
<input type="checkbox"/>	086F3	E0951 K0040 Adult Angle-Adj. w/Heel Loops ¹⁶ \$204.00

6. BACKREST / SEATING

BACKREST (16"-19", 8° w/Push Handles)

<input type="checkbox"/>	JAY J3 Back ¹⁵		
	086JB01	Please see the J3 Back order form for part #'s.	
<input type="checkbox"/>	086B12	Adjustable w/tall back canes ¹	\$15.00
<input type="checkbox"/>	086BU14	Omit Back Upholstery	-\$15.00

7. ARMRESTS

Flip Back Armrests

<input type="checkbox"/>		Fixed Height, Desk Length	Std
<input type="checkbox"/>	086AA5	Fixed Height, Full Length	NC
<input type="checkbox"/>	086AA6	E0973 Flip Back Height-Adjustable, Desk	\$317.00
<input type="checkbox"/>		E0973 Flip Back Height-Adjustable, Full	\$317.00

8. KEY OPTIONS AND ACCESSORIES

Backpack & Seat Pouch

<input type="checkbox"/>	086A1	Utility Back Pack (Black/Charcoal)	\$81.00
<input type="checkbox"/>	086A8	Seat Pouch (Black/Charcoal)	\$54.00

Wheelchair Tray Table - Tuffak

<input type="checkbox"/>	086A31	E0950 Medium (15"-17")	\$268.00
<input type="checkbox"/>	086A32	E0950 Large (18"-20")	\$268.00

Pelvic Positioning Belt

<input type="checkbox"/>	086A12	E0978 Aircraft Buckle (2" wide)	\$118.00
<input type="checkbox"/>	086A142	E0978 1.5" Auto Buckle	\$75.00

Miscellaneous

<input type="checkbox"/>	086A82	Caddy	\$54.00
<input type="checkbox"/>	086A51	E2207 Crutch Holder	\$209.00
<input type="checkbox"/>	086A83	Tool Kit	\$59.00

9. COLOR

Frame Color

<input type="checkbox"/>	086C6	Forest Green	Std
<input type="checkbox"/>	086C16	Silver Texture	NC
<input type="checkbox"/>	086C1	Blue	NC
<input type="checkbox"/>	086C2	Black	NC



The chart below specifies which casters and rear wheels are necessary to achieve the front seat heights listed below for **Standard** frames.

Seat Height	Rear Wheels				
	20" Low Poly	20" Full Poly/Pneu 22" Low Poly	22" Full Poly/Pneu	24" Low Poly	24" Full Poly/Pneu
15.0"	5" Ip	N/A	N/A	N/A	N/A
16.0"	5" Ip	5" Ip	N/A	N/A	N/A
16.5"	5" Ip		N/A	N/A	N/A
17.0"	5" Ip & 1-1/2" stem	5" Ip & 1-1/2" stem	5" Ip & 1-1/2" stem	5" Ip & 1-1/2" stem	N/A
17.5"	5" Ip & 1-1/2" stem or 8" Ip	5" Ip & 1-1/2" stem or 8" Ip	5" Ip & 1-1/2" stem or 8" Ip	5" Ip & 1-1/2" stem or 8" Ip	5" Ip & 1-1/2" stem or 8" Ip
18.0"	N/A	5" Ip & 1-1/2" stem	5" Ip & 1-1/2" stem	5" Ip & 1-1/2" stem	5" Ip & 1-1/2" stem
18.5"	N/A	8" Ip <small>*20" full poly not recommended</small>	8" Ip	8" Ip	8" Ip
19.0"	N/A	N/A	8" Ip & 1-1/2" stem or 8" pneu	8" Ip & 1-1/2" stem or 8" pneu	8" Ip & 1-1/2" stem or 8" pneu
19.5"	N/A	N/A	N/A	N/A	8" pneu & 1-1/2" stem
20.0"	N/A	N/A	N/A	N/A	8" Ip & 1-1/2" stem
20.5"	N/A	N/A	N/A	N/A	8" pneu & 1-1/2" stem

The chart below specifies which casters and rear wheels are necessary to achieve the front seat heights listed below for **Ultra Hemi** frames.

Seat Height	Rear Wheels				
	20" Low Poly	20" Full Poly/Pneu 22" Low Poly	22" Full Poly/Pneu	24" Low Poly	24" Full Poly/Pneu
13.0"	5" Ip	N/A	N/A	N/A	N/A
14.0"	5" Ip	5" Ip	N/A	N/A	N/A
14.5"	5" Ip & 1-1/2" stem	5" Ip & 1-1/2" stem	5" Ip & 1-1/2" stem	5" Ip & 1-1/2" stem	N/A
15.0"	8" Ip	8" Ip	8" Ip	8" Ip	N/A
15.5"	8" Ip or 5" Ip & 1-1/2" stem	8" Ip or 5" Ip & 1-1/2" stem	8" Ip or 5" Ip & 1-1/2" stem	8" Ip or 5" Ip & 1-1/2" stem	8" Ip or 5" Ip & 1-1/2" stem
16.0"	8" Ip or 8" pneu	8" Ip or 8" pneu	8" Ip or 8" pneu	8" Ip or 8" pneu	8" Ip or 8" pneu
16.5"	8" Ip & 1-1/2" stem	8" Ip & 1-1/2" stem	8" Ip & 1-1/2" stem	8" Ip & 1-1/2" stem	8" Ip & 1-1/2" stem
17.0"	8" Ip & 1-1/2" stem or 8" pneu	8" Ip & 1-1/2" stem or 8" pneu	8" Ip & 1-1/2" stem or 8" pneu	8" Ip & 1-1/2" stem or 8" pneu	8" Ip & 1-1/2" stem or 8" pneu
17.5"	N/A	1-1/2" stem & 8" Ip or 8" pneu	1-1/2" stem & 8" Ip or 8" pneu	1-1/2" stem & 8" Ip or 8" pneu	1-1/2" stem & 8" Ip or 8" pneu
18.0"	N/A	N/A	8" pneu & 1-1/2" stem	8" pneu & 1-1/2" stem	8" pneu & 1-1/2" stem