


## Arc Seat Pan



| Pan (Items 1 and 2) Part Numbers |                      |                |                 |                |                      |
|----------------------------------|----------------------|----------------|-----------------|----------------|----------------------|
|                                  | Frame Depth          | Frame Width    |                 |                |                      |
|                                  |                      | Small (10-12") | Medium (12-14") | Large (14-16") | Extra Large (16-18") |
| Item 1                           | Extra Short (12-15") | 113269         | 113270          | 113271         | 113272               |
|                                  | Short (14-17")       | 113273         | 113274          | 113275         | 113276               |
|                                  | Medium (16-19")      | 113277         | 113278          | 113279         | 113280               |
|                                  | Long (18-20")        | 113281         | 113282          | 113283         | 113284               |
| Item 2                           | All                  | 111430         | 111431          | 111432         | 111433               |

| Item Number | Part Number | Description                                            | Quantity |
|-------------|-------------|--------------------------------------------------------|----------|
| 1a, 2a, 3-8 | 113253      | Seat Pan Assembly, Extra Short Depth, 10-12" Wide, Arc | 1        |
| 1b, 2b, 3-8 | 113254      | Seat Pan Assembly, Extra Short Depth, 12-14" Wide, Arc | 1        |
| 1c, 2c, 3-8 | 113255      | Seat Pan Assembly, Extra Short Depth, 14-16" Wide, Arc | 1        |
| 1d, 2d, 3-8 | 113256      | Seat Pan Assembly, Extra Short Depth, 16-18" Wide, Arc | 1        |

| Item Number | Part Number | Description                                                       | Quantity |
|-------------|-------------|-------------------------------------------------------------------|----------|
| 1e, 2a, 3-8 | 113257      | Seat Pan Assembly, Short Depth, 10-12" Wide, Arc                  | 1        |
| 1f, 2b, 3-8 | 113258      | Seat Pan Assembly, Short Depth, 12-14" Wide, Arc                  | 1        |
| 1g, 2c, 3-8 | 113259      | Seat Pan Assembly, Short Depth, 14-16" Wide, Arc                  | 1        |
| 1h, 2d, 3-8 | 113260      | Seat Pan Assembly, Short Depth, 16-18" Wide, Arc                  | 1        |
| 1i, 2a, 3-8 | 113261      | Seat Pan Assembly, Medium Depth, 10-12" Wide, Arc                 | 1        |
| 1j, 2b, 3-8 | 113262      | Seat Pan Assembly, Medium Depth, 12-14" Wide, Arc                 | 1        |
| 1k, 2c, 3-8 | 113263      | Seat Pan Assembly, Medium Depth, 14-16" Wide, Arc                 | 1        |
| 1l, 2d, 3-8 | 113264      | Seat Pan Assembly, Medium Depth, 16-18" Wide, Arc                 | 1        |
| 1m, 2a, 3-8 | 113265      | Seat Pan Assembly, Long Depth, 10-12" Wide, Arc                   | 1        |
| 1n, 2b, 3-8 | 113266      | Seat Pan Assembly, Long Depth, 12-14" Wide, Arc                   | 1        |
| 1o, 2c, 3-8 | 113267      | Seat Pan Assembly, Long Depth, 14-16" Wide, Arc                   | 1        |
| 1p, 2d, 3-8 | 113268      | Seat Pan Assembly, Long Depth, 16-18" Wide, Arc                   | 1        |
| 1a, 4       | 113269      | Seat Pan with AB Tape, Front, Extra Short Depth, 10-12" Wide, Arc | 1        |
| 1b, 4       | 113270      | Seat Pan with AB Tape, Front, Extra Short Depth, 12-14" Wide, Arc | 1        |
| 1c, 4       | 113271      | Seat Pan with AB Tape, Front, Extra Short Depth, 14-16" Wide, Arc | 1        |
| 1d, 4       | 113272      | Seat Pan with AB Tape, Front, Extra Short Depth, 16-18" Wide, Arc | 1        |
| 1e, 4       | 113273      | Seat Pan with AB Tape, Front, Short Depth, 10-12" Wide, Arc       | 1        |
| 1f, 4       | 113274      | Seat Pan with AB Tape, Front, Short Depth, 12-14" Wide, Arc       | 1        |
| 1g, 4       | 113275      | Seat Pan with AB Tape, Front, Short Depth, 14-16" Wide, Arc       | 1        |
| 1h, 4       | 113276      | Seat Pan with AB Tape, Front, Short Depth, 16-18" Wide, Arc       | 1        |
| 1i, 4       | 113277      | Seat Pan with AB Tape, Front, Medium Depth, 10-12" Wide, Arc      | 1        |
| 1j, 4       | 113278      | Seat Pan with AB Tape, Front, Medium Depth, 12-14" Wide, Arc      | 1        |
| 1k, 4       | 113279      | Seat Pan with AB Tape, Front, Medium Depth, 14-16" Wide, Arc      | 1        |
| 1l, 4       | 113280      | Seat Pan with AB Tape, Front, Medium Depth, 16-18" Wide, Arc      | 1        |
| 1m, 4       | 113281      | Seat Pan with AB Tape, Front, Long Depth, 10-12" Wide, Arc        | 1        |
| 1n, 4       | 113282      | Seat Pan with AB Tape, Front, Long Depth, 12-14" Wide, Arc        | 1        |
| 1o, 4       | 113283      | Seat Pan with AB Tape, Front, Long Depth, 14-16" Wide, Arc        | 1        |

| Item Number | Part Number | Description                                                | Quantity |
|-------------|-------------|------------------------------------------------------------|----------|
| 1p, 4       | 113284      | Seat Pan with AB Tape, Front, Long Depth, 16-18" Wide, Arc | 1        |
| 2a          | 111430      | Seat Pan, Rear, 10-12 Wide                                 | 1        |
| 2b          | 111431      | Seat Pan, Rear, 12-14 Wide                                 | 1        |
| 2c          | 111432      | Seat Pan, Rear, 14-16 Wide                                 | 1        |
| 2d          | 111433      | Seat Pan, Rear, 16-18 Wide                                 | 1        |
| 3           | 111407      | Seat Pan Mount, Snap-on                                    | 4        |
| 4           | 001206      | AB Tape 3x3 Die Cut Tilt                                   | 2        |
| 5           | 102856      | Seat Pan Stop Collar                                       | 2        |
| 6           | 102634      | M6 X 1.0 X 16 BHCS BLK ZC                                  | 8        |
| 7           | 100746      | M6 Flat Washer Blk Zn                                      | 4        |
| 8           | 100658      | M6 Nylock Nut Blk Zn                                       | 8        |