


Account #:	_____	Purchase Order No.:	_____	Date:	_____
Contact:	_____	Mark For:	_____		
Bill To:	_____	Ship To:	_____		
Address:	_____	Address:	_____		
City/State/Zip:	_____	City/State/Zip:	_____		
Phone:	_____	Phone:	_____		
Email:	_____				

US

Effective: 01.10.2024 - Rev 3

Measurements



No.	Type	Measurement
9	Total height	
10	Shoulder height	
11	Chest height / Armpit height	
15	Top of pelvic height	
18	Chest circumference	
19	Pelvic circumference	

		0	1	2	3
	Width between handles	4 - 9¾"	7½ - 13¾"	9½ - 16½"	13 - 19"
	Height from floor to handles	11¾ - 24"	16½ - 28¾" 15 - 21¼"*	21¾ - 33"	29½ - 41¼"
	Max. user height	41¼"	45¼"	59"	70¾"
	Max. user mass	66 lbs	66 lbs	99 lbs	176 lbs
	Length	25¼"	28¼"	33"	34¼"
	Overall width	22¾"	24¾"	27¼"	27¼"
	Folded length	25½"	28½"	33¾"	39¼"
	Folded width	22¾"	24½"	27¼"	27¼"
	Folded height	11¾"	12¼"	10½"	11½"
	Minimum turning radius	19"	19¾"	22"	45"
	Mass device	13.4 lbs	15.7 lbs	15.7 lbs	26.7 lbs

*Mounted with the Starter Kit accessories. Further measurements at our website.

Features, incl. in standard price






- Width and height adjustable handles
- Easy to fold
- Lightweight durable aluminium frame – Sizes 0, 1 and 2
- Stable steel and aluminium frame – Size 3
- Upholstery and cover for top frame
- Solid front and rear wheels
- Swivel front wheels with directional stabilizer
- Tension adjustable rear wheels
- On/off anti-reverse wheel stop
- Padded Rear Stabilizer Bar w/ Tool Kit



Description	Measurements	Size	Article number	USD
-------------	--------------	------	----------------	-----

Main products

Crocodile Color and Size

	Crocodile, white Standard		Size 0	86800	718.00 <input type="checkbox"/>
	Crocodile, yellow Standard		Size 1	86801	792.00 <input type="checkbox"/>
	Crocodile, orange Standard		Size 2	86802	792.00 <input type="checkbox"/>
	Crocodile, metallic black Standard		Size 2	86852	792.00 <input type="checkbox"/>
	Crocodile, anthracite grey Standard		Size 3	86803	1,306.00 <input type="checkbox"/>



Description	Measurements	Size	Article number	USD
-------------	--------------	------	----------------	-----












Accessories








Mounting kits

	Crocodile Starter Kit <i>For size 1.</i>		One Size	86801-0	502.00 <input type="checkbox"/>
---	--	--	----------	---------	---------------------------------

Anti-tippers

	Anti-tip - Swing-away Anti-tip Swing-away, Pair, for size 0-2		Size 1	86815	171.00 <input type="checkbox"/>
	Anti-tip Swing-away, Pair, for size 3		Size 2	86815-3S	171.00 <input type="checkbox"/>
	Anti-tip - Fixed Anti-tip fixed, Pair - for size 1-2		Size 1	86815-F	82.00 <input type="checkbox"/>
	Anti-tip fixed, Pair - for size 3		Size 2	86815-3F	82.00 <input type="checkbox"/>

Description	Measurements	Size	Article number	USD
Arm supports				
	Forearm support w/handgrip			
	Forearm support, Pair - for size 0-1	L:5¾ x W:3"	Size 1	86858-1 462.00 <input type="checkbox"/>
	Forearm support, Pair - for size 2	L:7 x W:4"	Size 2	86859-1 462.00 <input type="checkbox"/>
	Forearm support, Pair - for size 3	L:7 x W:4"	Size 3	86860-1 462.00 <input type="checkbox"/>
	Padded forearm positioner strap*		One Size	
	Padded forearm positioner strap, EA			869370 73.00 <input type="checkbox"/>
	Padded forearm positioner strap, EA			869370 73.00 <input type="checkbox"/>
	<i>*For use with Forearm support w/ handgrip. For quantity 2 select both options above.</i>			
	Handle Knobs, Pair		One Size	86861 54.00 <input type="checkbox"/>
Brakes				
	Hand Brake Kit*			
	Hand Brake Kit, for size 1-2		Size 1	86834 256.00 <input type="checkbox"/>
	Hand Brake Kit, for size 3		Size 2	86873 256.00 <input type="checkbox"/>
	<i>*Hand brake eliminates the function of anti-reverse wheel stop</i>			
	Holder for Hand Brake*			
	Small, for size 1		Size 1	86874-1 359.00 <input type="checkbox"/>
	Medium, for size 2		Size 2	86874-2 405.00 <input type="checkbox"/>
	Large, for size 3		Size 3	86874-3 436.00 <input type="checkbox"/>
	<i>*Required with use of Hand brake and Forearm support. Allows for mounting of handbrake to Forearm support.</i>			
	Decelerating Hand Brake Kit*			
	Decelerating Hand Brake Kit, for sizes 1-2		Size 1	86834-DEC 256.00 <input type="checkbox"/>
	Decelerating Hand Brake Kit, for size 3		Size 2	86873-DEC 256.00 <input type="checkbox"/>
	<i>*Hand brake eliminates the function of anti-reverse wheel stop</i>			
	<i>*If ordering hand brake with anti-reverse wheel stop function, a custom order required</i>			
Back and pelvic supports				
	Surrounding Support*		One Size	868610 524.00 <input type="checkbox"/>
	<i>*Only available with size 0</i>			
	Pelvic Support w/ Laterals			
	Pelvic Support w/ Laterals, for sizes 0-2	3½" x 4¾" - 9"	Small	86808 281.00 <input type="checkbox"/>
	Pelvic Support w/ Laterals, for sizes 1-2	3½" x 7" - 11½"	Large	86809 281.00 <input type="checkbox"/>
	Lateral Strap*		One Size	905011 54.00 <input type="checkbox"/>
	<i>*For sizes 1 and 2. Back and side support size small and large</i>			
	Back support - Adjustable Padded*	9½" - 11½" x 7"		86886 541.00 <input type="checkbox"/>
	<i>*For Size 3</i>			
	Abdominal Support Belt*	13½" x 4"	-	85490-2 131.00 <input type="checkbox"/>
	<i>*For Back Support #86886, for size 3</i>			
	Cross bar for mounting back/pelvic supports, for Size 3		-	86880 104.00 <input type="checkbox"/>
	<i>*Required with all Back and Pelvic Supports for size 3.</i>			

Description	Measurements	Size	Article number	USD
Handles and push braces				
 Attendant Handle* *Not available with size 3.		One Size	86866-1	149.00 <input type="checkbox"/>
 Grip Bar Grip Bar, for size 0		Size 0	86879-0	138.00 <input type="checkbox"/>
Grip Bar, for size 1		Size 1	86879-1	159.00 <input type="checkbox"/>
Postural supports				
 Hip Supports, Pair* *Not available with size 3		One Size	86819	251.00 <input type="checkbox"/>
Seating				
Flip Down Seat White		Size 0	86820-01	182.00 <input type="checkbox"/>
Yellow*		Size 1	86820	182.00 <input type="checkbox"/>
Orange		Size 2	86821	182.00 <input type="checkbox"/>
Metallic Black		Size 2	86823	182.00 <input type="checkbox"/>
Anthracite Grey		Size 3	86863	182.00 <input type="checkbox"/>
*Not compatible with 86801-0, Crocodile Starter Kit				
 Lightweight Sling Seat* Lightweight Sling Seat, for size 0-1	23½" (Waist)	Size 1	86824-1	144.00 <input type="checkbox"/>
Lightweight Sling Seat, for size 2	27¼" (Waist)	Size 2	86825-1	144.00 <input type="checkbox"/>
*Requires Pelvic Support w/laterals 8680X and strap 905011 or 868610 for Size 0.				
 Sling seat Sling Seat, for size 1	23½" (Waist)	Size 1	86824	144.00 <input type="checkbox"/>
Sling Seat, for size 2	27¼" (Waist)	Size 2	86825	144.00 <input type="checkbox"/>
Bags and baskets				
 Bag Storage Bag, for sizes 0-2		One Size	86839	76.00 <input type="checkbox"/>
 Bag Storage Bag, for size 3		One Size	86843	225.00 <input type="checkbox"/>

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Prices in USD & is subject to change without notice

For Quotes and Orders please send completed Order Form to: sales.us@etac.com or fax: 1-800-760-0537

Etac
13137 Bleinheim Lane
Matthews, North Carolina 28105
+1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF5200-ENG-REV.03



Creating Possibilities