

Features



Padded armpad



Modern aesthetics with backrest detail



Easy to attach and remove legrests



Elevating Legrest Option

Essentials™ EX1

Technical data

For more information about this product, including the user manual, please visit our website: [www.invacareamerica.com](http://www.invacareamerica.com)

	Seat width	Overall width	Folded width	Seat depth	Overall depth (without riggings)	Overall height
Essentials™ EX1	20 inches (50.8 cm) 18 inches (45.7 cm) 16 inches (40.6 cm)	Seat width + 7.5 inches (19 cm)	11.8 inches (30 cm)	16 inches (40.6 cm)	31 inches (78.7 cm)	36 inches (91.4 cm)

	Seat to Floor Height	Back height	Rear wheel size	Weight	Shipping weight	Weight limit
Essentials™ EX1	19.5 inches (49.5 cm)	18 inches (46 cm)	24 inches (61 cm)	Chair only: 40 lbs (18.5 kg) Chair with front rigging: 44 lb (20 kg)	45.5 lb (20.6 kg)	250 lb (113.4 kg)

Item No.	Description
IEX166DASA	Essential EX1 16x16 Desk Arms Swing-Away Footrests
IEX186DASA	Essential EX1 18x16 Desk Arms Swing-Away Footrests
IEX106DASA	Essential EX1 20x16 Desk Arms Swing-Away Footrests
IEX166DAEL	Essential EX1 16x16 Desk Arms Swing-Away Legrests
IEX186DAEL	Essential EX1 18x16 Desk Arms Swing-Away Legrests
IEX106DAEL	Essential EX1 20x16 Desk Arms Swing-Away Legrests



Frame color

Invacare America  
One Invacare Way  
Elyria, Ohio 44035-4190  
(800) 333-6900

©2025 Invacare America. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare® America or its subsidiaries unless otherwise noted. Specifications are subject to change without notification. Form No. 24-036 / Job No. 250120 Rev 04.2025



Yes, you can.®



A durable and reliable

mobility solution




Yes, you can.®

# Introducing the new Essentials™ EX1

A lightweight mobility solution for increased independence.

Everyday mobility is made easier with the Essentials EX1. Designed for comfort, durability and transportation to help ensure you feel confident wherever you go.

Built with simplicity and functionality in mind, the Essentials EX1 is simple to go about your daily activities with smartly designed desk arms, swing-away legrests and a cross-folding frame.



## Reliable drive

The robust and durable frame is specifically designed to give you a reliable drive in most everyday environments to help achieve a more active lifestyle.



## Maximum Comfort

Benefiting from a subtly padded backrest and seat, the Essentials EX1 allows you to immerse into the chair so you feel comfortable and well supported. Also easy to maintain in a wipeable vinyl finish.



## Foldable & Compact

The lightweight frame and cross folding design helps to ensure folding is simple and effortless. Compact when folded also makes storage and transportation easy.

Black vinyl seat & back upholstery for easy cleaning

Desk length armpad gives easier access up to and under tables



Swing-away or elevating legrests\*



\*Product specifications may vary from those displayed.