



**ATTACHABLE PUSH BAR
WITH PADDED FOAM GRIP
MAKES IT EASY TO PUSH!**



**HANDY LEVER ON CHAIR
BACK MAKES IT EASY
AND CONVENIENT TO
FOLD BACK DOWN!**



ATC22-R

22" Bariatric Aluminum Transport Chair with 12" Rear "Flat Free" Wheels

ATC22-R Red 22" Width, Fold Down Back 1/cs

FEATURES:

- Lightweight - weighs only 33 lbs.
- Supports individuals up to 450 lbs.
- Attachable push bar with padded foam grip makes it easy to push.
- 12" rear flat-free wheels.
- Comes with swing-away footrests.
- Adjustable leg support.
- Comes with companion operated hand brakes and standard lever brakes.
- Fixed, full length arms.
- Padded armrests provide additional comfort.
- Heavy-duty, nylon reinforced upholstery with two carry handles.
- Comes with lever and companion activated wheel locks.
- Limited Lifetime Warranty.

**BACK FOLDS DOWN FOR
CONVENIENT STORAGE!**



SPECIFICATIONS

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/ Riggings	Product Weight	Carton Shipping Weight	Weight Limit
ATC22-R	30"	12"	37.5"	21"	16"	22"	18"	41"	33 lbs.	39 lbs.	450 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*

transport chairs

To order call toll free: 877.224.0946

drive™