


# (Discontinued) Ethos Transit



*Ethos Transit cannot be retrofitted. If needing to add transit, please contact Customer Service for a Replacement Seat Frame.*

Item Number	Part Number	Description	Quantity
1a	003846	M6 x 50 BHCS, Custom, Blk Zc <i>Used if Swing Away Arms are not present</i>	4
1b	003232	M6 x 55 BHCS Blk Zc <i>Used if Swing Away Arms are present</i>	4
2a	003848	Back Bracket Saddle 1 1/8" <i>Used if Swing Away Arms are not present</i>	2
3	003849	Saddle For Dimple 1 1/8"	4
4	109296	Back Bracket with Transit, Black, Rigid	2
6	100658	M6 Nylock Nut Blk Zn	4
7	100746	M6 Flat Washer Blk Zn	4

Item Number	Part Number	Description	Quantity
8	101853	M6 Washer Blk Zn	4