


# Rogue Transit



*This Rogue Transit is effective starting with serial number RG0016612. For all mounting hardware, please see Rogue Backrest Mount and Hardware page. Rogue Transit cannot be retrofitted. If needing to add transit, please see the Rogue Rigid Frame page or call Customer Service.*

Item Number	Part Number	Description	Quantity
2-6	112551	Kit Rogue Transit	1
2	003232	M6 x 55 BHCS Blk Zc	4
3	003715	Back Saddle 1 3/8" Tube	2
4	003716	Saddle 1 3/8" Tube	4
5	004434	Back Bracket with Transit, Black, Rigid	2
6	100658	M6 Nylock Nut Blk Zn	4