


(Discontinued 2) Ethos / Rogue Fixed Height Back Post with Non-Adjustable Height Rigidizer Bar



Effective on Rogue starting with serial number RG0016612. If retrofitting to this backrest on Rogue, a new backrest mount and hardware may be required. Tension Straps (#2) are based on seat width and used in combination with back upholstery.

Item Number	Part Number	Description	Quantity
1a	108722	Fixed Height Rigid Backrest 5", 12" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1b	108724	Fixed Height Rigid Backrest 6", 12" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1c	108726	Fixed Height Rigid Backrest 7", 12" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1

Item Number	Part Number	Description	Quantity
1d	108728	Fixed Height Rigid Backrest 8", 12" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1e	108730	Fixed Height Rigid Backrest 9", 12" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1f	108732	Fixed Height Rigid Backrest 10", 12" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1g	108734	Fixed Height Rigid Backrest 11", 12" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1h	108736	Fixed Height Rigid Backrest 12", 12" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1i	108737	Fixed Height Rigid Backrest 13", 12" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1j	108738	Fixed Height Rigid Backrest 14", 12" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1k	108739	Fixed Height Rigid Backrest 15", 12" Wide <i>Used with 17" back height with omit push handles</i>	1
1l	108740	Fixed Height Rigid Backrest 16", 12" Wide <i>Used with 18" back height with omit push handles</i>	1
1m	108741	Fixed Height Rigid Backrest 5", 13" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1n	108743	Fixed Height Rigid Backrest 6", 13" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1o	108745	Fixed Height Rigid Backrest 7", 13" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1p	108747	Fixed Height Rigid Backrest 8", 13" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1q	108749	Fixed Height Rigid Backrest 9", 13" Wide	1

Item Number	Part Number	Description	Quantity
		<i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	
1r	108751	Fixed Height Rigid Backrest 10", 13" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1s	108753	Fixed Height Rigid Backrest 11", 13" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1t	108755	Fixed Height Rigid Backrest 12", 13" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1u	108756	Fixed Height Rigid Backrest 13", 13" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1v	108757	Fixed Height Rigid Backrest 14", 13" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1w	108758	Fixed Height Rigid Backrest 15", 13" Wide <i>Used with 17" back height with omit push handles</i>	1
1x	108759	Fixed Height Rigid Backrest 16", 13" Wide <i>Used with 18" back height with omit push handles</i>	1
1y	108760	Fixed Height Rigid Backrest 5", 14" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1z	108762	Fixed Height Rigid Backrest 6", 14" Wide <i>Used with 10" back height with standard of low profile push handles or 12" back height with fold down push handles</i>	1
1aa	108764	Fixed Height Rigid Backrest 7", 14" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1ab	108766	Fixed Height Rigid Backrest 8", 14" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1ac	108768	Fixed Height Rigid Backrest 9", 14" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1ad	108770	Fixed Height Rigid Backrest 10", 14" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1

Item Number	Part Number	Description	Quantity
1ae	108772	Fixed Height Rigid Backrest 11", 14" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1af	108774	Fixed Height Rigid Backrest 12", 14" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1ag	108775	Fixed Height Rigid Backrest 13", 14" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1ah	108776	Fixed Height Rigid Backrest 14", 14" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1ai	108777	Fixed Height Rigid Backrest 15", 14" Wide <i>Used with 17" back height with omit push handles</i>	1
1aj	108778	Fixed Height Rigid Backrest 16", 14" Wide <i>Used with 18" back height with omit push handles</i>	1
1ak	108779	Fixed Height Rigid Backrest 5", 15" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1al	108781	Fixed Height Rigid Backrest 6", 15" Wide <i>Used with 10" back height with standard of low profile push handles or 12" back height with fold down push handles</i>	1
1am	108783	Fixed Height Rigid Backrest 7", 15" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1an	108785	Fixed Height Rigid Backrest 8", 15" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1ao	108787	Fixed Height Rigid Backrest 9", 15" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1ap	108789	Fixed Height Rigid Backrest 10", 15" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1aq	108791	Fixed Height Rigid Backrest 11", 15" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1ar	108793	Fixed Height Rigid Backrest 12", 15" Wide	1

Item Number	Part Number	Description	Quantity
		<i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	
1as	108794	Fixed Height Rigid Backrest 13", 15" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1at	108795	Fixed Height Rigid Backrest 14", 15" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1au	108796	Fixed Height Rigid Backrest 15", 15" Wide <i>Used with 17" back height with omit push handles</i>	1
1av	108797	Fixed Height Rigid Backrest 16", 15" Wide <i>Used with 18" back height with omit push handles</i>	1
1aw	108798	Fixed Height Rigid Backrest 5", 16" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1ax	108800	Fixed Height Rigid Backrest 6", 16" Wide <i>Used with 10" back height with standard of low profile push handles or 12" back height with fold down push handles</i>	1
1ay	108802	Fixed Height Rigid Backrest 7", 16" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1az	108804	Fixed Height Rigid Backrest 8", 16" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1ba	108806	Fixed Height Rigid Backrest 9", 16" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1bb	108808	Fixed Height Rigid Backrest 10", 16" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1bc	108810	Fixed Height Rigid Backrest 11", 16" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1bd	108812	Fixed Height Rigid Backrest 12", 16" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1be	108813	Fixed Height Rigid Backrest 13", 16" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1bf	108814	Fixed Height Rigid Backrest 14", 16" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1bg	108815	Fixed Height Rigid Backrest 15", 16" Wide <i>Used with 17" back height with omit push handles</i>	1
1bh	108816	Fixed Height Rigid Backrest 16", 16" Wide <i>Used with 18" back height with omit push handles</i>	1
1bi	108817	Fixed Height Rigid Backrest 5", 17" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1bj	108819	Fixed Height Rigid Backrest 6", 17" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1bk	108821	Fixed Height Rigid Backrest 7", 17" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1bl	108823	Fixed Height Rigid Backrest 8", 17" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1bm	108825	Fixed Height Rigid Backrest 9", 17" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1bn	108827	Fixed Height Rigid Backrest 10", 17" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1bo	108829	Fixed Height Rigid Backrest 11", 17" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1bp	108831	Fixed Height Rigid Backrest 12", 17" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1bq	108832	Fixed Height Rigid Backrest 13", 17" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1br	108833	Fixed Height Rigid Backrest 14", 17" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1bs	108834	Fixed Height Rigid Backrest 15", 17" Wide <i>Used with 17" back height with omit push handles</i>	1

Item Number	Part Number	Description	Quantity
1bt	108835	Fixed Height Rigid Backrest 16", 17" Wide <i>Used with 18" back height with omit push handles</i>	1
1bu	108836	Fixed Height Rigid Backrest 5", 18" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1bv	108838	Fixed Height Rigid Backrest 6", 18" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1bw	108840	Fixed Height Rigid Backrest 7", 18" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1bx	108842	Fixed Height Rigid Backrest 8", 18" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1by	108844	Fixed Height Rigid Backrest 9", 18" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1bz	108846	Fixed Height Rigid Backrest 10", 18" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1ca	108848	Fixed Height Rigid Backrest 11", 18" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1cb	108850	Fixed Height Rigid Backrest 12", 18" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1cc	108851	Fixed Height Rigid Backrest 13", 18" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1cd	108852	Fixed Height Rigid Backrest 14", 18" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1ce	108853	Fixed Height Rigid Backrest 15", 18" Wide <i>Used with 17" back height with omit push handles</i>	1
1cf	108854	Fixed Height Rigid Backrest 16", 18" Wide <i>Used with 18" back height with omit push handles</i>	1
1cg	108855	Fixed Height Rigid Backrest 5", 19" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1ch	108857	Fixed Height Rigid Backrest 6", 19" Wide <i>Used with 10" back height with standard of low profile push handles or 12" back height with fold down push handles</i>	1
1ci	108859	Fixed Height Rigid Backrest 7", 19" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1cj	108861	Fixed Height Rigid Backrest 8", 19" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1ck	108863	Fixed Height Rigid Backrest 9", 19" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1cl	108865	Fixed Height Rigid Backrest 10", 19" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1cm	108867	Fixed Height Rigid Backrest 11", 19" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1cn	108869	Fixed Height Rigid Backrest 12", 19" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1co	108870	Fixed Height Rigid Backrest 13", 19" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1cp	108871	Fixed Height Rigid Backrest 14", 19" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1cq	108872	Fixed Height Rigid Backrest 15", 19" Wide <i>Used with 17" back height with omit push handles</i>	1
1cr	108873	Fixed Height Rigid Backrest 16", 19" Wide <i>Used with 18" back height with omit push handles</i>	1
1cs	108874	Fixed Height Rigid Backrest 5", 20" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1ct	108876	Fixed Height Rigid Backrest 6", 20" Wide <i>Used with 10" back height with standard of low profile push handles or 12" back height with fold down push handles</i>	1
1cu	108878	Fixed Height Rigid Backrest 7", 20" Wide	1

Item Number	Part Number	Description	Quantity
		<i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	
1cv	108880	Fixed Height Rigid Backrest 8", 20" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1cw	108882	Fixed Height Rigid Backrest 9", 20" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1cx	108884	Fixed Height Rigid Backrest 10", 20" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1cy	108886	Fixed Height Rigid Backrest 11", 20" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1cz	108888	Fixed Height Rigid Backrest 12", 20" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1da	108889	Fixed Height Rigid Backrest 13", 20" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1db	108890	Fixed Height Rigid Backrest 14", 20" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1dc	108891	Fixed Height Rigid Backrest 15", 20" Wide <i>Used with 17" back height with omit push handles</i>	1
1dd	108892	Fixed Height Rigid Backrest 16", 20" Wide <i>Used with 18" back height with omit push handles</i>	1
2a	001749	Tension Strap Folding Handle 12-13 Wide	1
2b	001750	Tension Strap Folding Handle 14-15 Wide	1
2c	001751	Tension Strap Folding Handle 16-17 Wide	1
2d	001752	Tension Strap Folding Handle 18-19 Wide	1
2e	001753	Tension Strap Folding Handle 20-22 Wide	1
3	100804	Tube End Plug 1"	2
4	000073	M5 x 24 Threaded Barrel	2
5	100669	M5 x 10 BTN SHCS BLKZC w/Patch	4

Item Number	Part Number	Description	Quantity
6, 7a, 8	502602	Standard Push Handle Assembly, Short <i>Used with 9" - 18" back heights</i>	1
6, 7b, 8	107205	Low Profile Push Handle Assembly, Short <i>Used with 9" - 18" back heights</i>	1
7a	109432	Ergo Push Handle Grip 3/4" Standard with Instruction	1
7b	109430	Ergo Push Handle Grip 3/4" Low Profile with Instruction	1
8	100638	Tube End Plug 3/4"	1
9	000034	Sleeve 3/4" x 1 1/2"	2
10	100734	M5 x 16 PHMS ZC <i>Used to mount back upholstery and/or Fold Down Push Handle</i>	2
10, 11a, 12, 13	110529	Fold Down Push Handle Assembly, Fixed, 12" Back Height <i>Used with 12" back height</i>	1
10, 11b, 12, 13	110530	Fold Down Push Handle Assembly, Fixed, 13-18" Back Height <i>Used with 13" - 18" back heights</i>	1
12	001735	Folding Handle Assembly Rigid <i>This is to replace an existing Fold Down Push Handle. If wanting to retrofit to the Fold Down Push Handle, order one of the following complete assemblies based on the desired back height: 110529 or 110530</i>	1