

* Fab Angle (A) Frame Angle For use w/ Standard Sling								
Inset	0			1"			2	
Caster Setting	1	2	3	1	2	3	2	3
12" wide	114116S	114115S	125081S	114118S	114117S	125082S	N/A	N/A
13" wide	114128S	114127S	125087S	114130S	114129S	125088S	N/A	N/A
14" wide	114140S	114139S	125234S	114143S	114141S	125235S	114143S	125236S
15" wide	114155S	114154S	125243S	114157S	114156S	125244S	114158S	125245S
16" wide	114170S	114169S	125252S	114172S	114171S	125253S	114173S	125254S
17" wide	114185S	114184S	125261S	114187S	114186S	125262S	114188S	125263S
18" wide	114200S	114199S	125270S	114202S	114201S	125271S	114203S	125272S
19" wide	114215S	114214S	125138S	114217S	114216S	125139S	114218S	125140S
20" wide	114230S	114229S	125147S	114232S	114231S	125148S	114233S	125149S

* Fab Angle (B) Frame Angle For use w/ Standard Sling								
Inset	0			1"			2	
Caster Setting	1	2	3	1	2	3	2	3
12" wide	114120S	114119S	125083S	114122S	114121S	125084S	N/A	N/A
13" wide	114132S	114131S	125089S	114134S	114133S	125231S	N/A	N/A
14" wide	114145S	114144S	125237S	114147S	114146S	125238S	114148S	125239S
15" wide	114160S	114159S	125246S	114162S	114161S	125247S	114163S	125248S
16" wide	114175S	114174S	125255S	114177S	114176S	125256S	114178S	125257S
17" wide	114190S	114189S	125264S	114192S	114191S	125265S	114193S	125266S
18" wide	114205S	114204S	125273S	114207S	114206S	125274S	114208S	125134S
19" wide	114220S	114219S	125141S	114222S	114221S	125142S	114223S	125143S
20" wide	114235S	114234S	125150S	114237S	114236S	125151S	114238S	125152S

* Fab Angle (C) Frame Angle For use w/ Standard Sling								
Inset	0			1"			2	
Caster Setting	1	2	3	1	2	3	2	3
12" wide	114124S	114123S	125085S	114126S	114125S	125086S	N/A	N/A
13" wide	114136S	114135S	125232S	114138S	114137S	125233S	N/A	N/A
14" wide	114150S	114149S	125240S	114152S	114151S	125241S	114153S	125242S
15" wide	114165S	114164S	125249S	114167S	114166S	125250S	114168S	125251S
16" wide	114180S	114179S	125258S	114182S	114181S	125259S	114183S	125260S
17" wide	114195S	114194S	125267S	114197S	114196S	125268S	114198S	125269S
18" wide	114210S	114209S	125135S	114212S	114211S	125136S	114213S	125137S
19" wide	114225S	114224S	125144S	114227S	114226S	125145S	114228S	125146S
20" wide	114240S	114239S	125153S	114242S	114241S	125154S	114243S	125155S

When you ENTER the frame part number you will be prompted to provide the following SPECIFICATIONS:

- 1) Frame Color
- 2.) Frame Style
- 3.) Frame Depth
- 4.) Front Seat Height
- 5.) FS Dim
- 6.) DAB Option
- 7.) Q7 Wheel Lock Option
- 8.) Serial # w/ a -1 or -2 behind it depending on how many times the frame has been replaced.

*If making any OTHER CHANGES that will AFFECT STFH such as changing the caster size, fork size or rear wheel size; PLEASE go into QA-TEST (Kittens) reconfigure the chair with the new options. This is to verify that the new set up is compatible with the original chair and/or to capture all NO GO's that may exist. In MANY INSTANCES the FS Dimension or Frame Angle will NEED to be changed and seat to floor heights will increase or decrease, the dealer will need to be notified of the changes REQUIRED to get the new configuration.