

Raz-ART

Attendant Rotational Tilt

Mobile Shower Commode Chair



FEATURES

Patented Tilt Geometry:

· Dual arc, orbital tilt, minimizes center-of-mass movement throughout range

Extremely low tilt effort

• 40° tilt range

Small Footprint:

Provides unparalled manoeuvrability

Patent-Pending Tilt Lock:

Eliminates cable looping

Lightweight, invertible, aluminum frame:

Provides two tool-less height ranges totaling 7" (18 cm)



Flip-up arm supports with optional locks

Adjustable back angle

Compatible with all standardwidth Raz seats, which include the exclusive IPAS

Backrest with 3" (7.5 cm) fore/aft adjustment to vary seat depth

Built-in tilt angle indicator

Dual-locking (swivel and roll) 5" (125 mm) casters

Extensive range of options and accessories is available

ODEOUEIO ATIONIO	De ADT (A(Contest on the Contest De Contest TU)
SPECIFICATIONS	Raz-ART (Attendant propel Rotational Tilt)
WEIGHT CAPACITY	350 lbs (160 kg)
OVERALL WIDTH	22.5" (57 cm)
OVERALL LENGTH – with head support and standard foot supports in middle position	Upright: 35" (89 cm) / Tilted to 40°: 50.5" (128 cm)
SEAT WIDTH	18" (46 cm)
SEAT DEPTH – front of back post to front of seat	16" – 18" (40.5 – 46 cm) or 19" – 21" (48 – 53 cm) ¹
SEAT-TO-FLOOR HEIGHT – 1" (2.5 cm) increments – above standard molded seat	Lower Range: 20" – 23" (51 – 58.5 cm) Higher Range: 24" – 27" (61 – 68.5 cm)
CLEARANCE HEIGHT UNDER SEAT – 1" (2.5 cm) increments	
with pan brackets	Lower Range: 15" – 18" (38 – 46 cm) Higher Range: 19" – 22" (48 – 56 cm)
 without pan brackets 	Lower Range: 15.5" – 18.5" (41 – 49 cm) Higher Range: 19.5" – 22.5" (51 – 59 cm)
ARMREST HEIGHT – above standard molded seat	
– 18"H Fixed Back	9" (23 cm)
– 15"H Adjusta-Back	5.5" – 8.5" (14 – 21.5 cm) in 0.5" (1.25 cm) increments
– 17"H Adjusta-Back	6.5" – 9.5" (16.5 – 24 cm) in 0.5" (1.25 cm) increments
– 20"H Adjusta-Back	7.5" – 10.5" (19 – 27 cm) in 0.5" (1.25 cm) increments
– 24"H Adjusta-Back	8.5" – 11.5" (21.5 – 29 cm) in 0.5" (1.25 cm) increments
PRODUCT WEIGHT (with foot supports)	69 lbs (31 kg)
FRAME MATERIAL – base frame and seat frame FRAME MATERIAL – back frame and footrests HARDWARE MATERIAL	Powder-Coated Aluminum ² Stainless Steel ² Stainless Steel ²
¹ Ischial Pelvic Alignment System (IPAS) allows for the following s	seat adjustments: 2" (5 cm) of fore/aft, 1" (2.5cm) lateral, or 4.5° rotational

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.

MANUFACTURED BY: Raz Design Inc. 135 Railside Road, Toronto, Ontario Canada, M3A 1B7 T: 001.416.751.5678 E: info@razdesigninc.com

² Lifetime warranty against corrosion