


# Liberty Folding Base Cross Brace - Growth



Can you achieve your new width within the existing width range? If no, please choose the appropriate kit below. Kits below are preset to the widths shown.

Item Number	Part Number	Description	Quantity
1a, 2a, 3-23	110071	Cross Brace Assembly, Base, Small, 14W with Instruction <i>Used on 14 - 15" wide with folding frame. Has the ability to achieve all small width range settings (14" - 18").</i>	1
1a, 2a, 3-23	110072	Cross Brace Assembly, Base, Small, 16W with Instruction <i>Used on 16 - 18" wide with a small width range and folding frame. Has the ability to achieve all small width range settings (14" - 18").</i>	1
1b, 2b, 3-23	110073	Cross Brace Assembly, Base, Medium, 16W with Instruction <i>Used on 16 - 17" wide with medium width range and folding frame. Has the ability to achieve all medium width range settings (16" - 20").</i>	1
1b, 2b, 3-23	110074	Cross Brace Assembly, Base, Medium, 18W with Instruction	1

Item Number	Part Number	Description	Quantity
		<i>Used on 18 - 20" wide with medium width range and folding frame, Has the ability to achieve all medium width range settings (16" - 20").</i>	
1c, 2c, 3-23	110075	<b>Cross Brace Assembly, Base, Large, 18W with Instruction</b> <i>Used on 18 - 19" wide with large width range and folding frame. Has the ability to achieve all large width range settings (18" - 22").</i>	1
1c, 2c, 3-23	110076	<b>Cross Brace Assembly, Base, Large, 20W with Instruction</b> <i>Used on 20 - 22" wide with large width range and folding frame. Has the ability to achieve all large width range settings (18" - 22").</i>	1
1a, 2a, 3-8, 10, 11, 15, 16, 18-21	110077	<b>Cross Brace Assembly Non-Folding, Base, Small, 14W with Instruction</b> <i>Used on 14 - 15" wide with non-folding frame. Includes non-folding frame hardware. Has the ability to achieve all small width range settings (14" - 18").</i>	1
1a, 2a, 3-8, 10, 11, 15, 16, 18-21	110078	<b>Cross Brace Assembly Non-Folding, Base, Small, 16W with Instruction</b> <i>Used on 16 - 18" wide with a small width range and non-folding frame. Includes non-folding frame hardware. Has the ability to achieve all small width range settings (14" - 18").</i>	1
1b, 2b, 3-8, 10, 11, 15, 16, 18-21	110079	<b>Cross Brace Assembly Non-Folding, Base, Medium, 16W with Instruction</b> <i>Used on 16 - 17" wide with medium width range and non-folding frame. Includes non-folding frame hardware. Has the ability to achieve all medium width range settings (16" - 20"). Includes non-folding frame hardware.</i>	1
1b, 2b, 3-8, 10, 11, 15, 16, 18-21	110080	<b>Cross Brace Assembly Non-Folding, Base, Medium, 18W with Instruction</b> <i>Used on 18 - 20" wide with medium width range and non-folding frame. Includes non-folding frame hardware. Has the ability to achieve all medium width range settings (16" - 20"). Includes non-folding frame hardware.</i>	1
1c, 2c, 3-8, 10, 11, 15, 16, 18-21	110081	<b>Cross Brace Assembly Non-Folding, Base, Large, 18W with Instruction</b> <i>Used on 18 - 19" wide with large width range and non-folding frame. Includes non-folding frame hardware. Has the ability to achieve all large width range settings (18" - 22").</i>	1
1c, 2c, 3-8, 10, 11, 15, 16, 18-21	110082	<b>Cross Brace Assembly Non-Folding, Base, Large, 20W with Instruction</b> <i>Used on 20 - 22" wide with large width range and non-folding frame. Includes non-folding frame hardware. Has the ability to achieve all large width range settings (18" - 22").</i>	1
13, 22, 23	111130	<b>Folding Strap with Mount and Strap Slide, Liberty</b>	1