

## Pilot Plus

Cables Only			
Part Number:	139137		
Cable Length:	48" or 4'		
Part Number:	139138		
Cable Length:	61.5" or 5.125'		
Part Number:	139139		
Cable Length:	78.75" or 6.5'		

Unavailable / Obsolete Pilot				
Joystick Only				
Part Number: 126245				
Part Number:	126246			
Part Number:	118788			
Part Number:	118789			
Part Number:	139520			
Part Number:	139665			

VR2				
Joystick and Cable				
4 Button		6 Button		
Part Number:	105375	Part Number:	105376	
Plug:	Male	Plug:	Male	
Cable Length:	1200mm (47" or 4')	Cable Length:	1200mm (47" or 4')	
Obsolete VR2 - 4 Button		Obsolete VR2 - 6 Button		
Part Number:	101634	Part Number:	102174	
Part Number:	105248	Part Number:	105249	

## **R-Net**

CJSM 2 - BT and IR				
Joystick and Cable				
Part Number: 161664				
Plug: Male				
Cable Length:	1200mm (47" or 4')			
Obsolete CJSM 2				
Part Number:	131201			

When in doubt, order long male cable versions.

(Plug: Male, Cable Legth: 1200mm)

If the customer has a female joystick, a male joystick will replace the joystick and the extension cable.

CJSM 1 - Female Plug						
	All CJSM 1 have two Jack Ports on the bottom of the Joystick					
	Joystick and Cable					
STD		Lights w/ Pot and Toggle		Pot and Toggle		
Part Number:	108208	Part Number:	108213	Part Number:	108210	
Plug:	Female	Plug:	Female	Plug:	Female	
Cable Length:	111mm or 4"	Cable Length:	96mm or 3.75"	Cable Length:	96mm or 3.75"	
Obsolete CJSM 1		Obsolete CJSM 1 - Lights		Obsolete CJSM 1 - Toggle		
Part Number:	015571	Part Number:	108211	Part Number:	108209	
Part Number:	116525	Part Number:	108212			

CJSM 1 - Male Plug						
	All CJSM 1 have two Jack Ports on the bottom of the Joystick					
	Joystick and Cable					
STD		Lights		Pot and Toggle		
Part Number:	131263	Part Number:	131406	Part Number:	131397	
Plug:	Male	Plug:	Male	Plug:	Male	
Cable Length:	1200mm (47" or 4')	Cable Length:	1200mm (47" or 4')	Cable Length:	1200mm (47" or 4')	
Obsolete CJSM 1 - STD						
Part Number:	130392					
- Special Resistive Jack Input -						
If replacing this joystick and client has						
a 5-switch, must replace the 5-switch						
w/ a CTRL+5						

LED					
Joystick and Cable					
STD			Lights		
108207	Part Number:	131396	Part Number:	108206	
Female	Plug:	Male	Plug:	Female	
135mm or 5"	Cable Length:	1200mm (47" or 4')	Cable Length:	135mm or 5"	
Obsolete LED, STD, Female					
015570					
	108207 Female 135mm or 5" LED, STD, Female	STD  108207 Part Number: Female Plug: 135mm or 5" Cable Length: LED, STD, Female	Joystick and Cable  STD  108207 Part Number: 131396  Female Plug: Male  135mm or 5" Cable Length: 1200mm (47" or 4')  LED, STD, Female	Joystick and Cable  STD  108207 Part Number: 131396 Part Number:  Female Plug: Male Plug:  135mm or 5" Cable Length: 1200mm (47" or 4') Cable Length:  LED, STD, Female	