





The sweeping curves of the titanium dual-tube frame are the perfect fusion of striking design and driven performance. Lighter, more

efficient, and more elegant than any chair anywhere. As TiLite's flagship rigid chair, the TR is simply the best of the best.





TECHNICAL SPECIFICATIONS.

Frame Features

Frame Options

Transport Weight* Weight Limit Total Chair Width* Seat Width Front Seat Width (Seat Taper) Seat Depth Ergonomic Seat Front Seat Height **Rear Seat Height** Camber Front Frame Angle Footrest Width Options Center of Gravity Adjustment Front Wheels Rear Wheels Handrims Convenience Aesthetics

Titanium Dual-Tube Rigid TiFit Limited Adjustability **TR Black Edition** SuperLite Edition Armed Forces Edition Power Adaptable Frame Reinforced Frame 9.3 lbs. 265 lbs. 22" 12" - 20" 4 Options 12" - 20" 7 Options 15" - 21" 12.5" - 21" 0°, 2°, 4°, 6°, 8°, 12°, 15° or 18° 60°, 65°, 70°, 75°, 80°, 85° or 90° Standard Taper or V-Front End 5" 23 Options 67 Options 15 Options Folding Angle Adjustable Seat Back 6 Titanium Frame Finishes 19 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors Lifetime on Frame Limited Lifetime on Components

Warranty

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lbs.



TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit[®] chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

To learn more about TiLite options that can take your chair's performance to the next level, visit tilite.com/optimize.



TiLite | The Ultimate Ride™ 800.545.2266 | tilite.com | © 2014 TiLite, LLC. TiLite wheelchairs are shown with available options.