


# Rogue ALX Pelvic Positioning Belts - 2" Airline Buckle and Hardware (Formerly Tsunami)



Item Number	Part Number	Description	Quantity
1, 2, 3, 4	111403	Belt Pelvic Positioning 2" Airline Buckle 50", RG ALX/Rogue/LW XP	1
1	100594	Belt Pelvic Positioning 2" Airline Buckle 50"	1
2	101781	M6 x 50mm Partial Thread BHCS Blk ZC	2
3	100746	M6 Flat Washer Blk Zn	2
4	100658	M6 Nylock Nut Blk Zn	2