


Rogue XP Backrest Mount - Growth



Back Release Cable based on seat width The Back Release Cable length will measure 1" less than the specified seat width

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1a, 3 | 121180 | Back Release Cable 8" Wide Black with Lock Out Pins | 1 |
| 1b, 3 | 121179 | Back Release Cable 9" Wide Black with Lock Out Pins | 1 |
| 1c, 3 | 121181 | Back Release Cable 10" Wide Black with Lock Out Pins | 1 |
| 1d, 3 | 121182 | Back Release Cable, 11" Wide, Black with Lock Out Pins | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1e, 3 | 121183 | Back Release Cable, 12" Wide, Black with Lock Out Pins | 1 |
| 1f, 3 | 121184 | Back Release Cable, 13" Wide, Black with Lock Out Pins | 1 |
| 1g, 3 | 121185 | Back Release Cable, 14" Wide, Black with Lock Out Pins | 1 |
| 1h, 3 | 121186 | Back Release Cable, 15" Wide, Black with Lock Out Pins | 1 |
| 1i, 3 | 121187 | Back Release Cable, 16" Wide, Black with Lock Out Pins | 1 |
| 1j, 3 | 121188 | BACK RELEASE CABLE, 17" WIDE, BLACK with Lock Out Pins | 1 |
| 1k, 3 | 121189 | BACK RELEASE CABLE, 18" WIDE, BLACK with Lock Out Pins | 1 |
| 2 | 000965 | Backrest Cable Release Spring RG | 2 |
| 3 | 119660 | Lockout Pin | 2 |