


## Catalyst E Backrest - Angle Adjustable



| Item Number    | Part Number | Description                                    | Quantity |
|----------------|-------------|------------------------------------------------|----------|
| 1-3a, 4-10, 12 | 115834      | Angle Adjustable Backrest - Short (16")        | 1        |
| 1-3b, 4-10, 12 | 115835      | Angle Adjustable Backrest - Tall (20")         | 1        |
| 4 - 8, 10, 12  | 115836      | Angle Adjustable Backrest, Pivot Assembly      | 1        |
| 9, 11, 13, 14  | 115852      | Backrest Mount Hardware                        | 1        |
| 4              | 101593      | M6 x 25 BTN SHCS BLKZC                         | 3        |
| 5              | 111959      | Angle Adjustable Backrest Angle Plate Sleeve   | 2        |
| 6              | 111958      | Angle Adjustable Backrest Plate                | 1        |
| 7              | 112984      | Angle Adjustable Backrest Angle Plate Threaded | 1        |
| 8              | 100558      | M6 Nyl Ins Domed Cap L/N Zc                    | 1        |
| 9              | 100658      | M6 Nylock Nut Blk Zn                           | 2        |
| 10             | 111960      | Angle Adjustable Backrest Pivot Weldment       | 1        |
| 11             | 100860      | M6 X 35 Hhcs P/Thd Zc                          | 1        |
| 12             | 102625      | M6 X 30 BHCS BLK ZC                            | 1        |
| 13             | 100746      | M6 Flat Washer Blk Zn                          | 1        |
| 14             | 111965      | Frame Sleeve, Angle Adjustable Backrest        | 1        |