


Focus MAXX Spoke Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and Handrim kits will be listed. These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. For Focus chairs, new wheel locks might be needed when changing wheel size or style.

Item Number	Part Number	Description	Quantity
	109840	24" MAXX Spoke, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount	1
	109915	24" MAXX Spoke, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	119364	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	120112	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119365	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120113	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119366	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120114	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119367	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120115	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119368	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120116	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119369	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120117	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109841	24" MAXX Spoke, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	109842	24" MAXX Spoke, Full Poly Tire, Projection Handrim Kit, Standard Mount	1
	118199	24" MAXX Spoke, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount	1
	118215	24" MAXX Spoke, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118200	24" MAXX Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118201	24" MAXX Spoke, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount	1
	119479	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit	1

Item Number	Part Number	Description	Quantity
	120221	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119480	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120222	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119481	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120223	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119482	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120224	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119483	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120225	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119484	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120226	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	118112	24" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount	1
	118113	24" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	119470	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	120214	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119471	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120213	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119472	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120212	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1

Item Number	Part Number	Description	Quantity
	119467	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120211	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119468	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120210	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119469	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120209	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	118120	24" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118121	24" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount	1