


Rogue2 / Rogue XP Luggage Carrier



Used on Rogue XP chairs invoiced starting on 3/12/25

Item Number	Part Number	Description	Quantity
1 - 8	504465	Luggage Carrier Kit, Rogue/Rogue2/Rogue XP/RG ALX/LW XP	1
1	003430	Luggage Carrier Fork Assembly	2
2	003412	Luggage Carrier Clamp 1.375 <i>Used to mount Luggage Carrier to the frame</i>	4
3	003413	Luggage Carrier Barrel	2
4	100780	M6 x 25 SHCS ZC <i>Used with 003412</i>	2
5	100658	M6 Nylock Nut Blk Zn	2
6	100667	M6 x 12 BTN SHCS BLKZC w/Patch	4
7	003409	Luggage Carrier Clamp .75 <i>Used to mount Luggage Carrier to the footrest tubes/extension</i>	4
8	101845	M6 x 22 SHCS Blk Zn <i>Used with 003409</i>	2