


## Clik Camber Tube - Growth



Black Camber Tube kits based on seat width. Color Anodized Camber Tube kits used on all chair widths. New camber tube required when changing the degree of camber by more than 2 degrees. If Wheelchair Drive Ready Option is present, replacement parts need to be processed through Alber USA. If OAD and 5th Wheel are present, please select the appropriate OAD Camber Tube and Standard Camber Tube when growing chair in width.

| Item Number | Part Number | Description   | Quantity |
|-------------|-------------|---|----------|
| 1a          | 503238      | Camber Tube Rigid 8" Kit Black  | 1        |
| 1b          | 503239      | Camber Tube Rigid 9" Kit Black  | 1        |
| 1c          | 500720      | Camber Tube Rigid 10" Kit Black   | 1        |
| 1d          | 500721      | Camber Tube Rigid 11" Kit Black   | 1        |
| 1e          | 003369      | One Arm Drive Camber Tube 11" Kit, Clik / Rogue XP<br><i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i> | 1        |
| 1f          | 500722      | Camber Tube Rigid 12" Kit Black   | 1        |
| 1g          | 003370      | One Arm Drive Camber Tube 12" Kit, Clik / Rogue / Rogue XP  | 1        |

| Item Number | Part Number | Description   | Quantity |
|-------------|-------------|---|----------|
|             |             | <i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i>   |          |
| 1h          | 500723      | Camber Tube Rigid 13" Kit Black   | 1        |
| 1i          | 003371      | One Arm Drive Camber Tube 13" Kit, Klik / Rogue / Rogue XP<br><i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i> | 1        |
| 1j          | 500724      | Camber Tube Rigid 14" Kit Black   | 1        |
| 1k          | 003372      | One Arm Drive Camber Tube 14" Kit, Klik / Rogue / Rogue XP<br><i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i> | 1        |
| 1l          | 500725      | Camber Tube Rigid 15" Kit Black   | 1        |
| 1m          | 003373      | One Arm Drive Camber Tube 15" Kit, Klik / Rogue / Rogue XP<br><i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i> | 1        |
| 1n          | 500726      | Camber Tube Rigid 16" Kit Black   | 1        |
| 1o          | 003374      | One Arm Drive Camber Tube 16" Kit, Klik / Rogue / Rogue XP<br><i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i> | 1        |
| 1p          | 500727      | Camber Tube Rigid 17" Kit Black   | 1        |
| 1q          | 003708      | One Arm Drive Camber Tube 17" Kit, Klik / Rogue / Rogue XP<br><i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i> | 1        |
| 1r          | 500728      | Camber Tube Rigid 18" Kit Black   | 1        |
| 1s          | 003709      | One Arm Drive Camber Tube 18" Kit, Klik / Rogue / Rogue XP<br><i>Used if One Arm Drive is present. Does not include tube end plugs; add quantity 2 of 100532.</i> | 1        |
| 1t          | 501749      | Camber Tube Rigid Blue Anodized   | 1        |
| 1u          | 501750      | Camber Tube Rigid Bright Green Anodized   | 1        |
| 1v          | 501751      | Camber Tube Rigid Orange Anodized   | 1        |
| 1w          | 501752      | Camber Tube Rigid Pink Anodized   | 1        |
| 1x          | 501753      | Camber Tube Rigid Pewter Anodized   | 1        |
| 1y          | 501754      | Camber Tube Rigid Purple Anodized   | 1        |
| 1z          | 501755      | Camber Tube Rigid Red Anodized  | 1        |
| 1aa         | 501756      | Camber Tube Rigid Gold Anodized   | 1        |
| 2           | 100667      | M6 x 12 BTN SHCS BLKZC w/Patch  | 4        |

| Item Number | Part Number | Description             | Quantity |
|-------------|-------------|-------------------------|----------|
| 3           | 000066      | M6 x 27 Threaded Barrel | 2        |

# Color Anodizing Samples

