



ORDER FORM

**Q50 R carbon**

EIPW50

Customer	Ship To	
Account#	Name	
Date	Address	
PO#	City	
Buyer	State	Zip
Mark For	Phone	

Q50 R Carbon Totals	\$
Base Price	
Total of additional equipment	
Packaging and transport	
Total Price	

## Base Specifications

User Weight Capacity	300 lbs.
Weight (no batteries)	32 lbs.
Seat Size	Width 17", Depth 17"
Drive Wheel Position	Rear-Wheel Drive
Base Size	Width 22.4", Length 36.3"
Drive Control	nVR2 Brushless Controller
Driving Range	7.5 miles (15 miles with additional battery)
Drive Motors	Brushless motor / 3.7 mph
Back Height	18"
STF Height	19.5"
Drive Wheels	8" with flat-free inserts
Tires / Casters	7" Solid
Charger	2A Off-Board
Battery Size	10Ah Lithium Ion
Color	Carbon Fiber with 5 color accents
Ground Clearance	2.75"

## Additional Included Features and Accessories

Full Carbon Fiber
Quick-Release Folding
Lap Belt
Transport Strap
1.25" Thick Seat Cushion
Underseat Storage Bag
Storage for Additional Battery
Cup Holder
Off-Chair Battery Charging
Rear Running Lights
Depth Adjustable Joystick Mount
Flip-Up Footplate
Flip-Up Armrests
Breathable Back Upholstery

## OPTIONS

OPTION	HCPC	DESCRIPTION	RETAIL
<b>Model</b>			<b>\$</b>
500MW93	K0899	Q50 R Carbon	\$3,667.00
<b>Colors</b>			<b>\$</b>
500C200	-	Black	NC
500C201	-	Red	NC
500C202	-	Blue	NC
500C203	-	Yellow	NC
500C204	-	Purple	NC
500C205	-	Orange	NC
<b>Additional Battery</b>			<b>\$</b>
500AG43	-	Additional 10 Ah Lithium Ion Battery	\$607.00

Signature