


Focus CR Wheel Lock



Item Number	Part Number	Description	Quantity	Retail
1a, 2b, 3a, 4 - 10	501486	Wheel Lock Assembly Push Left Focus	1	\$48.66
1b, 2a, 3a, 4 - 10	501487	Wheel Lock Assembly Push Right Focus	1	\$48.66
1a, 1b, 2a, 2b, 3a, 4 - 10	501113	Wheel Lock Assembly Push PAIR Focus	1	\$97.32
2b, 3a, 4 - 10, 11a, 12-15	501488	Wheel Lock Assembly Push w/ Ext Left Focus	1	\$81.47
2a, 3a, 4 - 10, 11b, 12-15	501489	Wheel Lock Assembly Push w/ Ext Right Focus	1	\$81.47

Item Number	Part Number	Description	Quantity	Retail
2a, 2b, 3a, 4 - 10, 11a, 11b, 12-15	501114	Wheel Lock Assembly Push w/ Ext PAIR Focus	1	\$162.92
1c, 2b, 3b, 4 - 10	501490	Wheel Lock Assembly Pull Left Focus	1	\$48.66
1d, 2a, 3b, 4 - 10	501491	Wheel Lock Assembly Pull Right Focus	1	\$48.66
1c, 1d, 2a, 2b, 3b, 4 - 10	501115	Wheel Lock Assembly Pull PAIR Focus	1	\$97.32
2b, 3b, 4 - 10, 11c, 12-15	501492	Wheel Lock Assembly Pull w/ Ext Left Focus	1	\$81.47
2a, 3b, 4 - 10, 11d, 12-15	501493	Wheel Lock Assembly Pull w/ Ext Right Focus	1	\$81.47
2a, 2b, 3b, 4 - 10, 11c, 11d, 12-15	501116	Wheel Lock Assembly Pull w/ Ext PAIR Focus	1	\$162.92
2a	001138	Plate Wheel Lock Right Tilt	1	\$16.32
2b	001139	Plate Wheel Lock Left Tilt	1	\$16.32
3a	003769	Wheel Lock Mount Push Tilt <i>Replaced 001140</i>	1	\$16.81
3b	003770	Wheel Lock Mount Pull Tilt <i>Replaced 001238</i>	1	\$32.53
4	001141	Weld Nut NDM 06028 Zinc	1	\$2.63
5	101824	M6 x 30 HHCS Blk Zn	1	\$2.10
6	101823	M6 x 20 HHCS Blk Zn	1	\$1.05
7	100658	M6 Nylock Nut Blk Zn	2	\$1.05
8	101853	M6 Washer Blk Zn	1	\$0.53
9	100746	M6 Flat Washer Blk Zn	1	\$1.58
10	101835	M6 Lock Washer Blk Zn	1	\$0.53
12-15	100854	Brake Extension Replacement	1	\$12.19
13	101117	M5 x 22 BTN SHCS BLK ZC	1	\$1.58
14	101856	Flat Washer, M5, Blk Zn	1	\$0.53
15	100657	M5 Nylock Nut Blk Zn	1	\$0.53

Inverted axle plates are not used with wheels 20" or larger.

Wheel	Tire	Caster	Standard Axle Plate Seat Heights	Inverted Axle Plate Seat Heights
12" Mag	Pneumatic	5"	13.5" - 14.5"	17.5" - 19.5"
12" Mag	Pneumatic	6"	14.5"	17.5" - 19.5"
12" Mag	Pneumatic	6x2"	14.5"	17.5" - 19.5"
12" Mag	Pneumatic	7"	n/a	17.5" - 19.5"
12" Mag	Pneumatic	8"	n/a	17" - 19.5"
12" Mag	Pneumatic	8"x1.5"	n/a	17.5" - 19.5"
12" Mag	Pneumatic	8x2"	n/a	17.5" - 19.5"
12" Mag	Full Poly	5"	n/a	16.5" - 19"
12" Mag	Full Poly	6"	n/a	16.5" - 19"
12" Mag	Full Poly	6x2"	n/a	16.5" - 19"
12" Mag	Full Poly	7"	n/a	16.5" - 19"
12" Mag	Full Poly	8"	n/a	16." - 18.5"
12" Mag	Full Poly	8"x1.5"	n/a	16.5" - 19"
12" Mag	Full Poly	8x2"	n/a	17" - 18.5"
12" Mag w/ Drum Brake	Pneumatic	5"	13" - 14.5"	17.5" - 20"
12" Mag w/ Drum Brake	Pneumatic	6"	14" - 14.5"	17.5" - 19.5"
12" Mag w/ Drum Brake	Pneumatic	6x2"	14" - 14.5"	17.5" - 20"
12" Mag w/ Drum Brake	Pneumatic	7"	n/a	17.5" - 20"
12" Mag w/ Drum Brake	Pneumatic	8"	n/a	17.5" - 20"
12" Mag w/ Drum Brake	Pneumatic	8"x1.5"	n/a	17.5" - 19.5"
12" Mag w/ Drum Brake	Pneumatic	8x2"	n/a	17.5" - 19.5"
12" Mag w/ Drum Brake	Full Poly	5"	13" - 14.5"	17.5" - 19.5"
12" Mag w/ Drum Brake	Full Poly	6"	14" - 14.5"	17.5" - 19.5"
12" Mag w/ Drum Brake	Full Poly	6x2"	14" - 14.5"	17.5" - 19.5"
12" Mag w/ Drum Brake	Full Poly	7"	n/a	17.5" - 19.5"
12" Mag w/ Drum Brake	Full Poly	8"	n/a	17.5" - 19.5"
12" Mag w/ Drum Brake	Full Poly	8"x1.5"	n/a	17.5" - 19.5"
12" Mag w/ Drum Brake	Full Poly	8x2"	n/a	17.5" - 19.5"

Inverted axle plates are not used with wheels 20" or larger.

Wheel	Tire	Caster	Standard Axle Plate Seat Heights	Inverted Axle Plate Seat Heights
16" Mag	Pneumatic	5"	13" - 16"	19" - 20.5"
16" Mag	Pneumatic	6"	14" - 16"	19" - 20.5"
16" Mag	Pneumatic	6x2"	14" - 16"	19" - 20.5"
16" Mag	Pneumatic	7"	15" - 16"	19" - 20.5"
16" Mag	Pneumatic	8"	16"	19" - 20.5"
16" Mag	Pneumatic	8"x1.5"	15.5" - 16"	19" - 20.5"
16" Mag	Pneumatic	8x2"	16"	19" - 20.5"
16" Mag	Full Poly	5"	13" - 16.5"	19" - 20.5"
16" Mag	Full Poly	6"	14" - 16.5"	19" - 20.5"
16" Mag	Full Poly	6x2"	14" - 16.5"	19" - 20.5"
16" Mag	Full Poly	7"	15" - 16.5"	19.5" - 20.5"
16" Mag	Full Poly	8"	16"	19" - 20.5"
16" Mag	Full Poly	8"x1.5"	15.5" - 16.5"	19" - 20.5"
16" Mag	Full Poly	8x2"	16" - 16.5"	19" - 20.5"
16" Mag	Low Poly	5"	13" - 15"	18" - 20"
16" Mag	Low Poly	6"	14" - 15"	18" - 20"
16" Mag	Low Poly	6x2"	14" - 15"	18" - 20.5"
16" Mag	Low Poly	7"	n/a	18" - 20.5"
16" Mag	Low Poly	8"	n/a	18" - 20.5"
16" Mag	Low Poly	8"x1.5"	n/a	18" - 20"
16" Mag	Low Poly	8x2"	n/a	18" - 20"
16" Mag w/ Drum Brake	All	5"	13.5" - 16"	19" - 20.5"
16" Mag w/ Drum Brake	All	6"	14" - 16"	18.5" - 20.5"
16" Mag w/ Drum Brake	All	6x2"	14" - 16"	18.5" - 20.5"
16" Mag w/ Drum Brake	All	7"	15" - 15.5"	18.5" - 20.5"
16" Mag w/ Drum Brake	All	8"	n/a	18.5" - 20.5"
16" Mag w/ Drum Brake	All	8"x1.5"	15.5" - 16"	18.5" - 20.5"
16" Mag w/ Drum Brake	All	8x2"	16"	18.5" - 20.5"