


# Focus CR Fixed Height with Adjustable Handle Backrest - Growth



Has your width changed? If yes, verify the current growth strut will work or select appropriately below. If chair was invoiced prior to 12/18/23, an upper back tube assembly is required. If chair was invoiced after to 12/18/23, upper back tube assembly parts can now be purchased individually

Item Number	Part Number	Description	Quantity
1a, 2	503107	Growth Strut Assembly Short, Focus/Flip	1

Item Number	Part Number	Description	Quantity
		<i>Used on chairs 14"-17" wide</i>	
1b, 2	503108	Growth Strut Assembly Long, Flip/Focus <i>Used on chairs 18"-22" wide</i>	1
1c, 2	114289	Growth Strut Assembly Extra Long, Focus <i>Used on chairs 23"-24" wide</i>	1