



Seat-2-GO Positioning Seat

PORTABLE • COMFORTABLE • CONVENIENT

Take it with you everywhere you go!!

A lightweight soft padded positioning seat for children that need support while sitting on the floor, at the table or at play. It is simple to use in a home, school or when traveling

CRS 2000 Small, 1/cs
CRS 3000 Medium, 1/cs

FEATURES:

- Contoured seat prevents child from sliding forward and improves seating posture.
- 95° seat to back angle.
- Double strap trunk support and belt.
- Fastening straps secure seat to any chair and can be stored in the back pocket when not in use.
- Can be used as an insert for a single or double stroller.
- 1 year limited warranty.



**WEIGHS ONLY
5 LBS!**



positioning
seating

Wenzel^{ite} Re/hab

To order call toll free:
New York: 877-224-0946 • California: 800-371-2266

drive
MEDICAL DESIGN & MANUFACTURING

Seat2GO

Positioning Seat

ACCESSORIES:

Headrest

CRS 8000 Headrest for CRS 2000 & CRS 3000, 1/bx

Height adjustable headrest.

Abductor

CRS 8025 Abductor for CRS 2000 & CRS 3000, 1/bx

Depth adjustable abductor.



positioning

seating

Wenzel^{ite} Rehab

To order call toll free:
New York: 877.224.0946 • California: 800.371.2266

ITEM #	WIDTH	DEPTH	HEIGHT	WEIGHT CAPACITY	PRODUCT WEIGHT	SHIP WEIGHT
CRS 2000	8"	9"	12"	50 lbs.	5 lbs.	8 lbs.
CRS 3000	12"	12"	14"	75 lbs.	7 lbs.	11 lbs.

drive
MEDICAL DESIGN & MANUFACTURING



WZ 33