


MAXX Mag Wheel / Tire Assemblies



These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or receivers may be needed. For the Catalyst 4 or Liberty, new axle plates may be needed. For CR45 if changing wheel size and Attendant Foot Lock or T-Arms are present, spacer 112503 is required which can be found on the CR45 Axle Plate page See chair specific axles page for axle and receiver utilization.

Item Number	Part Number	Description	Quantity
1a	116870	22" MAXX Mag Wheel with Pneumatic w/ Airless Insert Tire, Black	1
1b	116869	24" Wheel MAXX Mag 5-Spoke, Pneumatic w/ Airless Insert Tire, Black	1