



Rogue2 Footrest - Tubular with Carbon Fiber Cover




12" Rigid	10.1"
13" Rigid	11.1"
14" Rigid	12.1"
15" Rigid	13.1"
16" Rigid	14.1"
17" Rigid	15.1"
18" Rigid	16.1"

Not available with V-Taper options. Sizes based on straight frame. If changing from a high mount footrest to a frame inserted footrest, please see the "Rogue2 Rigid Frame" page for necessary hardware.

Item Number	Part Number	Description	Quantity
1a, 2a	125158	Tubular Footrest w/ Carbon Fiber Footplate Cover Assembly 12" Rigid	1
1b, 2a	125159	Tubular Footrest w/ Carbon Fiber Footplate Cover Assembly 13" Rigid	1
1c, 2a	125153	Tubular Footrest w/ Carbon Fiber Footplate Cover Assembly 14" Rigid	1
1d, 2b	125154	Tubular Footrest w/ Carbon Fiber Footplate Cover Assembly 15" Rigid	1
1e, 2b	125155	Tubular Footrest w/ Carbon Fiber Footplate Cover Assembly 16" Rigid	1
1f, 2c	125156	Tubular Footrest w/ Carbon Fiber Footplate Cover Assembly 17" Rigid	1
1g, 2c	125157	Tubular Footrest w/ Carbon Fiber Footplate Cover Assembly 18" Rigid	1
3, 4	501739	Footrest Protector w/ Instructions Kit Left Rigid	1
3, 4	501740	Footrest Protector w/ Instructions Kit Right Rigid	1
3, 4	501741	Footrest Protector w/ Instructions Kit PAIR Rigid	1
1a	000214	Footrest Tubular 12" Rigid	1
1b	000215	Footrest Tubular 13" Rigid	1
1c	000216	Footrest Tubular 14" Rigid	1
1d	000217	Footrest Tubular 15" Rigid	1
1e	000218	Footrest Tubular 16" Rigid	1
1f	000219	Footrest Tubular 17" Rigid	1
1g	000220	Footrest Tubular 18" Rigid	1

Distance between Point A and Point B
B = Footplate Length

NOTE: The yardstick on the top plane of the seat frame must lie flat against frame. The yardstick going from the footplate to the second yardstick must be square against the hanger tube.



Distance between Point A and Point B
B = Footplate Length

NOTE: The yardstick on the top plane of the seat frame must lie flat against frame. The yardstick going from the footplate to the second yardstick must be square against the front tube.



Distance between Point A and Point B
B = Footplate Length

NOTE: The yardstick on the top plane of the hanger must lie flat against hanger. The yardstick going from the footplate to the second yardstick must be square against the hanger tube.

