



Liberty Transit Wheels - Folding Frame

Rear Seat Height	Item 1 Left Side	Item 1 Right Side	Item 2
14"	113951	113950	111298
14.5"	113951	113950	111299
15" - 15.5"	113953	113952	111300
16" - 16.5"	113955	113954	111300
17" - 18"	113955	113954	111301



Not available with Long Seat Depth Range

Item Number	Part Number	Description	Quantity
1a, 1b, 2a, 3 - 5	115588	Liberty Transit Wheels Kit, Small/Extra Small (PAIR), Short/Medium Depth Range, Folding <i>Used with 14" rear seat height</i>	1
1a, 1b, 2b, 3 - 5	115589	Liberty Transit Wheels Kit, Small/Small (PAIR), Short/Medium Depth Range, Folding <i>Used with 14.5" rear seat height</i>	1
1c, 1d, 2c, 3 - 5	115590	Liberty Transit Wheels Kit, Medium/Medium (PAIR), Short/Medium Depth Range, Folding <i>Used with 15" - 15.5" rear seat heights</i>	1
1e, 1f, 2c, 3 - 5	115591	Liberty Transit Wheels Kit, Large/Medium (PAIR), Short/Medium Depth Range, Folding <i>Used with 16" - 16.5" rear seat heights</i>	1
1e, 1f, 2d, 3 - 5	115592	Liberty Transit Wheels Kit, Large/Large (PAIR), Short/Medium Depth Range, Folding <i>Used with 17" - 18" rear seat heights</i>	1
1a, 2a, 3 - 5	115593	Liberty Transit Wheel Kit, Small/Extra Small, Left Side, Short/Medium Depth Range, Folding <i>Used with 14" rear seat height</i>	1
1a, 2b, 3 - 5	115594	Liberty Transit Wheel Kit, Small/Small, Left Side, Short/Medium Depth Range, Folding <i>Used with 14.5" rear seat height</i>	1
1c, 2c, 3 - 5	115595	Liberty Transit Wheel Kit, Medium/Medium, Left Side, Short/Medium Depth Range, Folding <i>Used with 15" - 15.5" rear seat heights</i>	1

Item Number	Part Number	Description	Quantity
1e, 2c, 3 - 5	115596	Liberty Transit Wheel Kit, Large/Medium, Left Side, Short/Medium Depth Range, Folding <i>Used with 16" - 16.5" rear seat heights</i>	1
1e, 2d, 3 - 5	115597	Liberty Transit Wheel Kit, Large/Large, Left Side, Short/Medium Depth Range, Folding <i>Used with 17" - 18" rear seat heights</i>	1
1b, 2a, 3 - 5	115598	Liberty Transit Wheel Kit, Small/Extra Small, Right Side, Short/Medium Depth Range, Folding <i>Used with 14" rear seat height</i>	1
1b, 2b, 3 - 5	115599	Liberty Transit Wheel Kit, Small/Small, Right Side, Short/Medium Depth Range, Folding <i>Used with 14.5" rear seat height</i>	1
1d, 2c, 3 - 5	115600	Liberty Transit Wheel Kit, Medium/Medium, Right Side, Short/Medium Depth Range, Folding <i>Used with 15" - 15.5" rear seat heights</i>	1
1f, 2c, 3 - 5	115601	Liberty Transit Wheel Kit, Large/Medium, Right Side, Short/Medium Depth Range, Folding <i>Used with 16" - 16.5" rear seat heights</i>	1
1f, 2d, 3 - 5	115602	Liberty Transit Wheel Kit, Large/Large, Right Side, Short/Medium Depth Range, Folding <i>Used with 17" - 18" rear seat heights</i>	1
1a	113951	Liberty Transit Wheel Upper Assembly, Small, Left <i>Used with 14" - 14.5" rear seat heights</i>	1
1b	113950	Liberty Transit Wheel Upper Assembly, Small, Right <i>Used with 14" - 14.5" rear seat heights</i>	1
1c	113953	Liberty Transit Wheel Upper Assembly, Medium, Left <i>Used with 15" - 15.5" rear seat heights</i>	1
1d	113952	Liberty Transit Wheel Upper Assembly, Medium, Right <i>Used with 15" - 15.5" rear seat heights</i>	1
1e	113955	Liberty Transit Wheel Upper Assembly, Large, Left <i>Used with 16" - 18" rear sea heights</i>	1
1f	113954	Liberty Transit Wheel Upper Assembly, Large, Right <i>Used with 16" - 18" rear sea heights</i>	1
2a	111298	Liberty Transit Wheel Bottom Assembly, Extra Small, Short/Medium Depth Range <i>Used with 14" rear seat height</i>	1
2b	111299	Liberty Transit Wheel Bottom Assembly, Small, Short/Medium Depth Range <i>Used with 14.5" rear seat height</i>	1
2c	111300	Liberty Transit Wheel Bottom Assembly, Medium, Short/Medium Depth Range <i>Used with 15" - 16.5" rear seat heights</i>	1
2d	111301	Liberty Transit Wheel Bottom Assembly, Large, Short/Medium Depth Range <i>Used with 17" - 18" rear seat heights</i>	1
3	100658	M6 Nylock Nut Blk Zn	1

Item Number	Part Number	Description	Quantity
4	102072	Washer Nylon .234 ID x .406 OD x .062 Thick	1
5	100704	M6 x 90 CLS8.8 DIN931 HHCS ZC	1

Liberty FT Transit Wheels	
Initial Rear Seat Height	Height Adjustment Range
14"	N / A
14.5"	N / A
15"	15" - 15.5"
15.5"	15.5" - 16.5"
16"	15.5" - 16.5"
16.5"	16.5" - 18"
17"	16.5" - 18"
17.5"	17.5" - 18"
18"	17.5" - 18"

See Order Form for No-Go's