



Raz-AT600

MSRP – USA – USD
Effective Aug 28, 2023

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-AT600 Heavy Duty Mobile Shower Commode Chair with stainless steel frame and components, 25° tilt range, Fixed Back, four dual-locking 6" (150mm) casters, flip-up arm supports, heavy duty foot supports, commode pan, head support with small pad and the proprietary Weight Assist Spring System (WASS). Also included is a 20"H Multi-Position Back with 2" fore/aft adjustment range, two back angles (90°/97°) and tension-adjustable vinyl upholstery. **User weight capacity 605 lbs.** For help with justification, please see our [Evaluation and Medical Justification Form](#).

Need help with this order form? Call Raz customer service at (877) 720-5678

1. Chair Codes by size – **Must Choose One**

Z36022	22"W Seat and Back Frames	\$5,995
Z36024	24"W Seat and Back Frames	\$6,295
Z36026	26"W Seat and Back Frames	\$6,595
Z36028	28"W Seat and Back Frames	\$6,995
Z36030	30"W Seat and Back Frames – Requires custom seat	\$7,595

2. User Weight Range

Z36WR450	WASS User Weight Range 350-450 lbs.	Std	Z36WR600	WASS User Weight Range 400-605 lbs.	\$250
----------	-------------------------------------	-----	----------	-------------------------------------	-------

3. Seats – **Must Choose One**

Molded and Visco Foam Interface Seats – Seat Support Kit (Z4SSKAT) required if seat is narrower than seat frame.

ZMS23FRO	Molded 23"W x 18"D – Front Access	Fits 22"/24"W ¹ frames only	N/C
ZMS26FRO	Molded 26"W x 18"D – Front Access	Fits 26"/28"W ¹ frames only	N/C
ZVFI23FRO	Visco Foam Interface 23"W x 18"D – Front Access	Fits 22"/24"W ¹ frames only	\$625
ZVFI26FRO	Visco Foam Interface 26"W x 18"D – Front Access	Fits 26"/28"W ¹ frames only	\$745
ZSSKT	Seat Support Kit for Raz-AT600 – for seats that are 2" narrower than seat frames		\$145

¹ Required for seats that are narrower than seat frames

Custom Seats

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. [See Custom Seat Assessment Form](#) – Resources tab at razdesigninc.com

ZCSEATHD	Custom Dual-Visco Foam Seat – HD (enter specs below) Select if both width and depth are ≤ 28"	\$1,095
ZCSEATXHD	Custom Dual-Visco Foam Seat – XHD" (enter specs below) Select if width and/or depth is ≥ 28"	\$1,295

Base Frame Width: ___ 22" / 24" / 26" / 28" / 30"
Custom Seat Width: ___ 22" / 24" / 26" / 28" / 30" / 32" (cannot be more than 2" wider than seat frame width)
Seat Depth: ___ 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" (≥ 21" requires custom longer frame)
Access: ___ FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged
Aperture Width: ___ 6.5" recommended

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom seat widths can be up to 2" wider than seat the frame width to provide additional thigh support without increasing overall width of the chair.
- The Ischial Pelvic Alignment System (IPAS) will allow 2" fore/aft adjustment of the seat position (and therefore seat depth and aperture position) on all seats.
- Access for Rear-Open seats is limited due to rear seat support tube.

4. Symphony Back (Recessed solid back that replaces the standard Tension-Adjustable Vinyl Upholstery [recommended for infection-control environments]). All Symphony Back panels are 13"H.

ZSB22	Symphony Back – 22"W	\$465
ZSB24	Symphony Back – 24"W	\$485
ZSB26	Symphony Back – 26"W	\$505
ZSB28	Symphony Back – 28"W	\$575
ZSB30	Symphony Back – 30"W	\$675

5. Arm Support Accessories

ZALKHD	Arm Support Locks – HD (pr)	\$95	ZPAML	Pivoting Arm Mount – Left (for ZFAPLL / ZTRL)	\$165
ZFAPLL	Long Molded Flat Armpad ¹ – Left (14"L x 3.25"W)	\$90	ZPAMR	Pivoting Arm Mount – Right (for ZFAPLR / ZTRR)	\$165
ZFAPLR	Long Molded Flat Armpad ¹ – Right (14"L x 3.25"W)	\$90	ZAPSB_ _ *	Anterior Posture Support Bar * Enter frame width to complete code (22" / 24")	\$245
ZTRL	Molded Arm Trough ¹ – Left (18.5"L x 4.5"W)	\$110	ZAPSB_ _ *	Anterior Posture Support Bar * Enter frame width to complete code (26" / 28" / 30")	\$275
ZTRR	Molded Arm Trough ¹ – Right (18.5"L x 4.5"W)	\$110			
¹ May require Pivoting Arm Mount (ZPAM) for arm support to flip up					

6. Foot and Leg Supports

ZFSHD	Adjustable Footrests HD (pr) – 16"-20.5"	Std	ZELRL	Elevating Legrest ¹ – Left (incl ZADML)	\$495
ZFET2	Footrest Extension Tube (pr) – 2" longer	\$50	ZELRR	Elevating Legrest ¹ – Right (incl ZADMR)	\$495
ZADM	Angle/Depth Adj Footplates (pr)	\$120	ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 2.5" to accommodate deeper seats and reduced knee flexion range. Specify 1 or 2 _	\$95
ZFET2A	Footrest Extension Tube (pr) for ZADM – 2" longer	\$50	ZRLSL	Universal Residual Limb Support – Left	\$195
ZFBT	Flip-Back Foot Support (1-piece) with calf pads – 16"-20.5" length adjustment	\$455	ZRLSR	Universal Residual Limb Support – Right	\$195
			¹ ZELR compatible with adjusted seat depths ≤ 20"D		

7. Foot and Leg Support Accessories

ZFSLK	Foot Support Lock (pr)	\$95	ZNC	Neoprene Covers for Standard Footplates (pr)	\$90
ZCS__*	Snap-On Calf Strap – Bodypoint Aeromesh * Enter frame width to complete code (22" / 24")	\$105	ZNCA	Neoprene Covers for Angle/Depth Adjustable Footplates	\$90
ZCS__*	Snap-On Calf Strap – Bodypoint Aeromesh * Enter frame width to complete code (26" / 28" / 30")	\$245	ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$90
ZCSIC__*	Snap-On Infection Control Calf Strap * Enter frame width to complete code (22" / 24" / 26" / 28" / 30")	\$125	ZFOR	Forward Offset Foot Support Receiver - Required when seats extend beyond front edge of frame [provides 2.5" offset] Specify 1 or 2 ____	\$95
ZHS2224	H-Strap – for 22"/24"W frames – hook and loop attachment	\$110	ZLMOR	Lateral / Medial Offset Foot Support Receiver – moves foot support 2" laterally or medially. Specify 1 or 2 ____	\$95
ZAHL	Adjustable Heel Loops (pr)	\$60			

8. Positioning Pads – Pad Sizes: Small 5.5" x 3.5" / Large – 6.75" x 4.3".

Head Supports					
ZHSP2	Head Support Pad – Large (12" x 6")				\$50
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment				\$75
Lateral Thoracic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For up to 4.5" medial offset, specify ZLAEP extension plate.					
ZLTS2L	Lateral Thoracic Support w/ Large Pad – Left				\$225
ZLTS2R	Lateral Thoracic Support w/ Large Pad – Right				\$225
Leg Adductors					
ZLGT2L	Leg Adductor w/ Large Pad – Left				\$225
ZLGT2R	Leg Adductor w/ Large Pad – Right				\$225
Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors					
These plates allow for an additional 1.5" medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors					
ZLAEP	Lateral/Adductor Extension Plate (ea)				\$28
Lateral Thoracic Support	Left	Right			
Leg Adductor	Left	Right			
Hip Adductor	Left	Right			
Total Quantity ____					
Medial Thigh Supports (Pommel)					
ZMTST2	Medial Thigh Support w/ Large Pads				\$495

9. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts		
ZPBLHDT	Pelvic Belt – Large (seat frame widths 22"-24", 2-piece w/ side-release buckle)	\$65	ZCBL	Chest Belt – Large ² (back frame width 22"-24", 1-piece w/ hook and loop closure)	\$45
ZPBXLHDT	Pelvic Belt – XLarge (seat frame widths 24"-26", 2-piece w/ side-release buckle)	\$75	ZCBXL	Chest Belt – XLarge ² (back frame widths 24"-26", 1-piece w/ hook and loop closure)	\$55
ZPBXXLHDT	Pelvic Belt – XXL Large (seat frame widths 28"-30", 2-piece w/ side-release buckle)	\$110	ZCBXXL	Chest Belt – XXL Large ² (back frame widths 26"-30", 1-piece w/ hook and loop closure)	\$65
ZPBICLHDT	Infection Control Pelvic Belt – Large ¹ (seat frame widths 22"-24", 2-piece w/ side-release buckle)	\$115	ZCBICL	Infection Control Chest Belt – Large ^{1,2} (back frame width 22"-24", 1-piece w/ side-release buckle)	\$85
ZPBICXLHDT	Infection Control Pelvic Belt – XLarge ¹ (seat frame widths 24"-26", 2-piece w/ side-release buckle)	\$125	ZCBICXL	Infection Control Chest Belt – XLarge ^{1,2} (back frame widths 24"-26", 1-piece w/ side-release buckle)	\$95
ZPBICXXLHDT	Infection Control Pelvic Belt – XXL Large ¹ (seat frame widths 26"-30", 2-piece w/ side-release buckle)	\$135	ZCBICXXL	Infection Control Chest Belt – XLarge ^{1,2} (back frame widths 26"-30", 1-piece w/ side-release buckle)	\$105
ZPBBPLHDT	Pelvic Belt – Bodypoint Aeromesh – Large	\$170	ZCBBPL	Chest Belt – Bodypoint Aeromesh – Large	\$245
ZPBBPXLHDT	Pelvic Belt – Bodypoint Aeromesh – XLarge	\$235	ZCBBPXL	Chest Belt – Bodypoint Aeromesh – XLarge	\$195
Recommended sizes for Bodypoint belts	Pelvic belt	ZPBBPLHDT	22"-26"W base frame	ZPBBPXLHDT	26"-30"W base frame
	Chest belt	ZCBBPL	22"-26"W back frame	ZCBBPXL	26"-30"W back frame
¹ Non-absorbent belts – recommended for infection-control environments. ² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCB). When the chair is supplied with a Symphony Back, Chest Belt loops are provided.					

10. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$67
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$125
ZSC23HD	Splash Curtain ¹ for 23"W HD Molded and VFI Seats (surrounds entire aperture, raises seat 1/2").	\$345
ZSC26HD	Splash Curtain ¹ for 26"W HD Molded and VFI Seats (surrounds entire aperture, raises seat 1/2").	\$375
ZSCCHD	Custom Splash Curtain for Custom Seats – Raises seat 1/2"	\$365
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$345
ZBUDBH	Bard urine drainage bag holder for 22" / 24" / 26" / 28" / 30"W frames	\$155
¹ Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for pricing.		

11. Other

ZDCLT	Directional Caster Lock – fits on left rear caster (facilitates steering)	\$165	ZPSK1	Commode Pan Spacer Kit ² – 1" lower	\$200
ZATRR	Anti-Tippers – Rear (pr)	\$345	ZPSK2	Commode Pan Spacer Kit ² – 2" lower	\$215
ZATTFR	Anti-Tippers – Front (pr)	\$345	ZPSK3	Commode Pan Spacer Kit ² – 3" lower	\$230
ZTLS	Tilt Limit Spacers – 1/2" spacers snap onto gas struts to limit tilt (1/2" reduces tilt range approx. 5°. Can be stacked for additional tilt limiting. Quantity ____	N/C	ZPSK4	Commode Pan Spacer Kit ² – 4" lower	\$245
ZSPKTHD24	AT600 Self-Propel Kit – 24" wheels – includes rear anti-tippers	\$1,425	¹ Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Increases height required to clear toilet.		

12. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	–\$20	Z911	Omit Molded Head Support Pad	–\$20
Z907HD	Omit HD Foot Supports – if no other Foot / Leg Support Option is selected	–\$45	Z912	Omit Complete Head Support Assembly	–\$50

Some assembly of Raz products is required.

13. Notes:**Specifications****Raz-AT600 (Attendant Tilt)**

Weight Capacity	605 lbs
Seat Depth (with Standard Molded Seat)	18"-20"
Overall Width	Seat frame width plus 4.5" (with Standard Arm Support Pads)
Toilet Clearance Width	17"
Overall Length (with foot supports and headrest adjusted to middle positions)	
- at 0° tilt and standard foot supports	42"
- at 25°	52"
Clearance Height (under seat)	with 6" casters
- with pan brackets	17.75"-18.75"
- without pan brackets	18.5"-19.5"
Arm Support Height - above molded seat	9.25"
Seat-to-Floor Height (with molded seat)	22.25"
Back Height (Molded Seat to top of upholstery)	20"
Product Weight (with 22"W seat frame)	84 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)