


## Rogue2 22" MAXX Mag Wheel / Tire / Handrim Kits



*Not all components of the Wheel, Tire and Handrim kits will be listed. These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or receivers may be needed. For the Catalyst 4, new axle plates may be needed. See Rogue2 Axles page for axle and receiver utilization.*

Item Number	Part Number	Description	Quantity
	118218	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118218_ _ _	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	118219	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118219_ _ _	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	118220	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118220_ _ _	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	118221	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118221_ _ _	22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	118222	22" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118222_ _ _	22" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount, Colored	1
	118223	22" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118223_ _ _	22" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount, Colored	1
	118224	22" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount	1
	118224_ _ _	22" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount, Colored	1
	118225	22" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount	1
	118225_ _ _	22" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount, Colored	1
	119527	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119527_ _ _	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119528	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119528_ _ _	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119529	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119529_ _ _	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119531	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1

Item Number	Part Number	Description	Quantity
	119531_ _ _	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119532	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	119532_ _ _	22" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1
	118122	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118122_ _ _	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	118123	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118123_ _ _	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	118124	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118124_ _ _	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	118125	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118125_ _ _	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	118126	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118126_ _ _	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount, Colored	1
	118127	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118127_ _ _	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount, Colored	1
	118128	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount	1
	118128_ _ _	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount, Colored	1
	118129	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Close Mount	1
	118129_ _ _	22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Close Mount, Colored	1
	119521	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit	1

Item Number	Part Number	Description	Quantity
	119521_ _ _	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119522	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119522_ _ _	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119523	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119523_ _ _	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119525	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119525_ _ _	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119526	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	119526_ _ _	22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1
	117682	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117682_ _ _	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117683	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117683_ _ _	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117684	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117684_ _ _	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117685	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117685_ _ _	22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117686	22" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117686_ _ _	22" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount, Colored	1
	117687	22" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount	1
	117687_ _ _	22" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount, Colored	1
	117688	22" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	117688_ _ _	22" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Standard Mount, Colored	1
	117689	22" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Close Mount	1
	117689_ _ _	22" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Close Mount, Colored	1
	119497	22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119497_ _ _	22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119498	22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119498_ _ _	22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119499	22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119499_ _ _	22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119501	22" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119501_ _ _	22" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119502	22" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	119502_ _ _	22" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1
	117720	22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117720_ _ _	22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117721	22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117721_ _ _	22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117722	22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117722_ _ _	22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117723	22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117723_ _ _	22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117724	22" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117724_ _ _	22" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Standard Mount, Colored	1
	117725	22" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Close Mount	1
	117725_ _ _	22" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Close Mount, Colored	1

Item Number	Part Number	Description	Quantity
	119503	22" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119503_ _ _	22" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119504	22" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119504_ _ _	22" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119505	22" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119505_ _ _	22" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119507	22" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119507_ _ _	22" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119508	22" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	119508_ _ _	22" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1
	117754	22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117754_ _ _	22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount, Colored	1
	117755	22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117755_ _ _	22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount, Colored	1
	117756	22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117756_ _ _	22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip, Colored	1
	117757	22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117757_ _ _	22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip, Colored	1
	117758	22" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117758_ _ _	22" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount, Colored	1
	117759	22" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Close Mount	1
	117759_ _ _	22" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Close Mount, Colored	1
	119509	22" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119509_ _ _	22" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119510	22" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1

Item Number	Part Number	Description	Quantity
	119510_ _ _	22" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119511	22" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119511_ _ _	22" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119513	22" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119513_ _ _	22" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119514	22" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	119514_ _ _	22" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1
	117800	22" MAXX Mag, Mountain Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117801	22" MAXX Mag, Mountain Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117802	22" MAXX Mag, Mountain Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117803	22" MAXX Mag, Mountain Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117804	22" MAXX Mag, Mountain Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117805	22" MAXX Mag, Mountain Tire, Plastic Coated Handrim Kit, Close Mount	1
	117806	22" MAXX Mag, Mountain Tire, Projection Handrim Kit, Standard Mount	1
	117807	22" MAXX Mag, Mountain Tire, Projection Handrim Kit, Close Mount	1
	119515	22" MAXX Mag Mountain Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	119515_ _ _	22" MAXX Mag Mountain Tire Ergonomic Standard Handrim Omit Thumb Grip Kit Colored	1
	119516	22" MAXX Mag Mountain Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	119516_ _ _	22" MAXX Mag Mountain Tire Ergonomic Standard Handrim Regular Thumb Grip Kit Colored	1
	119517	22" MAXX Mag Mountain Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	119517_ _ _	22" MAXX Mag Mountain Tire Ergonomic Standard Handrim Extra Thumb Grip Kit Colored	1
	119519	22" MAXX Mag Mountain Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	119519_ _ _	22" MAXX Mag Mountain Tire Ergonomic Lite Handrim Regular Thumb Grip Kit Colored	1
	119520	22" MAXX Mag Mountain Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	119520_ _ _	22" MAXX Mag Mountain Tire Ergonomic Lite Handrim Extra Thumb Grip Kit Colored	1