


Flip for x:panda Transit



Item Number	Part Number	Description	Quantity
1-10	121066	Transit Backrest Hardware Kit, Single Side, Flip for x:panda	1
1-14	121067	Transit Retrofit Kit, Flip for x:panda	1
2	100746	M6 Flat Washer Blk Zn	2

Item Number	Part Number	Description	Quantity
3	001239	P-Clip 3/16"	2
4	100658	M6 Nylock Nut Blk Zn	8
5	112491	Saddle, Dual, 1" x 2.0SP	4
6	002738	Spacer .313 OD x .243 ID x 1.68 long	2
7	001897	Washer, Dual, PSA	4
9	101781	M6 x 50mm Partial Thread BHCS Blk ZC	6
10	003232	M6 x 55 BHCS Blk Zc	2
12	101733	M6 x 35 CLS10.9 DIN7991 FSHCS BLKZC	2
13	002086	Spacer, .38ODx.24IDx.31LG	2
14	100558	M6 Nyl Ins Domed Cap L/N Zc	2