



285

Available in
Black or White

Rehab Shower/Commode Chair-Tilt

- Packages Now Available
- Stainless Steel Frame
- Several New Seat Options Available
- Tilt Adjustment
- Four Frame Sizes:
16"wx18"D, 18"wx18"D, 20"wx18"D, and 22"wx18"D
- Seven Seat Sizes:
16"wx18"D, 18"wx18"D, 18"wx20"D, 20"wx18"D, 20"wx20"D,
22"wx18"D, and 22"wx20"D

The 285 is a gravity assist tilt-in-space shower/commode chair. Available in multiple frame widths from 16"-22" and the option of 4 casters or 24" rear wheels make this a versatile chair for the home or rehab environment.

Multiple features combine to provide support, comfort, and pressure relief along with the other benefits of tilt-in-space technology.



FRAMES



P13635
16"Wx18"D
Stainless Steel
Frame



P13636
18"Wx18"D
Stainless Steel
Frame



P13637
20"Wx18"D
Stainless Steel
Frame



P13638
22"Wx18"D
Stainless Steel
Frame



60800-2
5" Dual Locking
Casters
(pair)



P61096
24" Wheels
(pair)



P60549
24" Pneumatic
Street Tread
Wheels
(pair)

AXLES

WHEEL LOCKS

FRONT CASTERS



P61059
24" Flat Free Street
Tread Wheels
(pair)



P61104
24" PVC Handrims
(replaces standard
handrims, pair)



P13688
Bolted Axles with
Tipping Levers
(pair)



P60532
Wheel Locks
(pair)



P13135
Wheel Locks with
Extended Handle
(pair)



60841-2
5" Casters
(pair)



60800-2
5" Dual Locking
Casters
(pair)

SEATS



P61124
16"Wx18"D
Ensolite® 1" foam
Front/Rear Open Seat
A-4.75"xB-7.75" opening



P60993
18"Wx18"D
Comfortuff 1.5" waterfall foam
Front/Rear Open Seat
A-5"xB-8" opening



P60974
18"Wx18"D
Comfortuff sloped foam
Front Open Seat
A-4"xB-6.5" opening



P60802
18"Wx18"D
Ensolite® 1.25" foam
Front Open Ring Seat
A-4.5"xB-7.75" opening



P61086
18"Wx18"D
Ensolite® 2" foam
Solid Seat



P13771
18"Wx20"D
Ensolite® 1" foam
Front Open Seat
A-4.5"xB-7.5" opening



P61038
20"Wx18"D
Comfortuff
1.5" waterfall foam
Front Open Seat
A-5"xB-8" opening



P60903
20"Wx18"D
Ensolite® 1" foam
Front/Rear Open Seat
A-4.5"xB-7.75" opening



PP61122
20"Wx18"D
Ensolite® 1" foam
Solid Seat




P13772
20"Wx20"D
Ensolite® 1.5" foam
Race Track Oval Seat
A-7.75"xB-12.75" opening



P13773
20"Wx20"D
Ensolite® 1" waterfall foam
Front/Rear Open Seat
A-4.5"xB-7.75" opening



P60905
22"Wx18"D
Ensolite® 1" foam
Front/Rear Open Seat
A-4.5"xB-7.75" opening



P13774
22"Wx20"D
Ensolite® 1.5" foam
Race Track Oval Seat
A-7.75"xB-12.75" opening



P13775
22"Wx20"D
Ensolite® 1" waterfall foam
Front/Rear Open Seat
A-4.5"xB-7.75" opening

ARMRESTS



P13611
Left Hand Armrest
with Pad



P13612
Right Hand Armrest
with Pad



P13613
Left Hand Armrest
with Arm Trough



P13614
Right Hand Armrest
with Arm Trough



P13615
Left Hand
Height-Adjustable
Armrest



P13616
Right Hand
Height-Adjustable
Armrest

LEG RESTS

FOOT PLATES



PL13581
Left Gooseneck
Leg Rest



PR13581
Right Gooseneck
Leg Rest



P13582
Left Elevating
Leg Rest



P13583
Right Elevating
Leg Rest



PL11099
Left Foot Plate



PR11099
Right Foot Plate



PL11100
Left Foot Plate with
Heel Loop

FOOT PLATES

BACKS



PR11100
Right Foot Plate
with Heel Loop



P61137
Left Oversized
Aluminum
Foot Plate



P61138
Right Oversized
Aluminum
Foot Plate



PL11101
Left Padded
Foot Plate



PR11101
Right Padded
Foot Plate



P61139
16"W Adjustable
Sling Back



P11106
16"W Solid Back
Full Length w/
Center Pad



P11110
16"W Solid Back



P60877
18"W Adjustable
Sling Back



P11107
18"W Solid Back



P60878
20"W Adjustable
Sling Back



P11108
20"W Solid Back



P61140
22"W Adjustable
Sling Back



P11109
22"W Solid Back

BELTS



PW60884
Waist Belt
(48" circumference with
12" velcro overlap)



PT60910
Trunk Belt
(48" circumference with
12" velcro overlap)



P60879
Calf Strap



P11093
Bodypoint®
Trunk Belt



P11102
Bodypoint®
Trunk Belt (pair)



P11116
Bodypoint® Trunk
Belt & Waist Belt



P11112
Bodypoint® XL
Trunk Belt

POSITIONING AIDS




P11113
Belt w/Quick
Release Knobs



P61081
Bodypoint®
Calf Strap



P11114
Bodypoint®
Knee Belt



PL13574
Left Hand
Swing-Away
Adjustable Lateral



PR13574
Right Hand
Swing-Away
Adjustable Lateral



PL13618
Left Hand
Arm-Mounted
Hip Guide



PR13618
Right Hand
Arm-Mounted
Hip Guide

PANS

HEAD SUPPORTS



P61006
Abductor/Deflector
Small



P13224
Pan/Hanger for
Ring and
Contoured Seats



P13423
Pan/Hanger for
Front/Rear, & Race
Track Oval Seats




P11105
Head Support
Small Foam
(only available for
P11104)




P13621
Head Support
Large

Front Open Seat
Opening Measurements



Race Track Oval
Opening Measurements



full option details at
activeaid.com



285 with Support Options



PKG207 Package #5 - 18"W x 18"D

SPECIFICATIONS

Weight Capacity	350 lbs.
Maximum Tilt	30°
Top of Seat to Top of Back	22"
Floor to Top of Seat	25"
Seat Width	16", 18", 20" and 22"
Seat Depth	18" and 20"
Floor to Bottom of Seat Frame (without pan)	19.75"
Width Between Arm Pads	20.25"
Top of Seat to Top of Armrest Pads	7.5"
Overall Width	24.5"
Maximum Back Slack	3.5"
Headrest Height from Seat to Bottom of Cervical Support	20.25" - 28.25"
Minimum Overall length with Footrest – No Tilt	47"
Maximum Overall Length with Footrest – Full Tilt – Headrest in	61"
Legrest Range	15.5"-21.5"
Under Seat Width Clearance	17.5"