


# Liberty Drum Brake - Growth



*Did your depth change? If yes, evaluate if a different size cable is necessary.*

Item Number	Part Number	Description	Quantity
1a	112928	Cable, Tilt w/Drum, Long, Liberty	2
1b	112929	Cable, Tilt w/Drum, XL	2
1c	112930	Cable, Tilt w/Drum, XXL	2
1d	112931	Cable, Tilt w/Drum, XXXL	2
1e	112568	Cable, Tilt w/Drum, 4XL, Liberty	2
1f	112569	Cable, Tilt w/Drum, 5XL, Liberty	2

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
2a	002471	Cable, Drum Brake, XL	2
2b	002472	Cable, Drum Brake, XXL	2
2c	111547	Cable, Drum Brake, 3XL	2
2d	111548	Cable, Drum Brake, 4XL	2
2e	111549	Cable, Drum Brake, 5XL	2

Liberty FT Drum Brake Cable Usage							
Backrest Type & Height Range	Seat Depth						
	14"	15"	16"	17"	18"	19"	20"
Height Adjustable Straight, Short Height Adjustable 8° Bend, Short	XL (002471)	XL (002471)	2XL (002472)	2XL (002472)	3XL (111547)	3XL (111547)	4XL (111548)
Height Adjustable Straight, Medium Height Adjustable 8° Bend, Medium Stroller Handle, Short	XL (002471)	2XL (002472)	2XL (002472)	3XL (111547)	3XL (111547)	4XL (111548)	4XL (111548)
Height Adjustable Straight, Tall Height Adjustable 8° Bend, Tall	2XL (002472)	2XL (002472)	3XL (111547)	3XL (111547)	4XL (111548)	4XL (111548)	5XL (111549)
Stroller Handle, Tall	2XL (002472)	3XL (111547)	3XL (111547)	4XL (111548)	4XL (111548)	5XL (111549)	5XL (111549)

Liberty FT Tilt Cable Usage			
Backrest Type & Height Range	Seat Depth Range		
	Short (14" - 16")	Medium (16" - 18")	Long (18" - 20")
Height Adjustable Straight, Short Height Adjustable 8° Bend, Short	Long (112928)	XL (112929)	XXL (112930)
Height Adjustable Straight, Medium Height Adjustable 8° Bend, Medium Stroller Handle, Short	XL (112929)	XXL (112930)	XXXL (112931)
Height Adjustable Straight, Tall Height Adjustable 8° Bend, Tall	XXL (112930)	XXXL (112931)	4XL (112568)
Stroller Handle, Tall	XXXL (112931)	4XL (112568)	5XL (112569)