

## Jay Lite Cushion Measuring Guide

NOTE: Measurements are taken from the bottom of the foam base and have a +/- 0.25 tolerance.

Part Number	Actual Width	Actual Depth
JL1414	13.25"	14.50"
JL1416	13.25"	16.50"
JL1418	13.25"	18.50"
JL1515	14.25"	15.50"
JL1516	14.25"	16.50"
JL1517	14.25"	17.50"
JL1518	14.25"	18.50"
JL1615	15.25"	15.50"
JL1616	15.25"	16.50"
JL1617	15.25"	17.50"
JL1618	15.25"	18.50"
JL1620	15.25"	20.50"
JL1715	16.25"	15.50"
JL1716	16.25"	16.50"
JL1717	16.25"	17.50"
JL1718	16.25"	18.50"
JL1816	17.25"	16.50"
JL1817	17.25"	17.50"
JL1818	17.25"	18.50"
JL1820	17.25"	20.50"
JL1916	18.25"	16.50"
JL1917	18.25"	17.50"
JL1918	18.25"	18.50"
JL1919	18.25"	19.50"
JL1920	18.25"	20.50"
JL2016	19.25"	16.50"
JL2018	19.25"	18.50"
JL2020	19.25"	20.50"