


# One Piece Flip Up Footplate



Assemblies are based on seat width. If Tilt chair has Offset Frame, add 2" to seat width Only compatible with extension tube hangers.

| Item Number  | Part Number | Description                                       | Quantity |
|--------------|-------------|---|----------|
| 3-19         | 112894      | Hardware Assembly for One Piece Flip Up Footplate | 1        |
| 1a, 2a, 3-19 | 500517      | Footrest Assembly One Piece Flip 14" Wide         | 1        |
| 1a, 2b, 3-19 | 500518      | Footrest Assembly One Piece Flip 15" Wide         | 1        |
| 1a, 2c, 3-19 | 500519      | Footrest Assembly One Piece Flip 16" Wide         | 1        |

| Item Number  | Part Number | Description                               | Quantity |
|--------------|-------------|---|----------|
| 1b, 2d, 3-19 | 500520      | Footrest Assembly One Piece Flip 17" Wide | 1        |
| 1b, 2e, 3-19 | 500521      | Footrest Assembly One Piece Flip 18" Wide | 1        |
| 1b, 2f, 3-19 | 500522      | Footrest Assembly One Piece Flip 19" Wide | 1        |
| 1b, 2g, 3-19 | 500523      | Footrest Assembly One Piece Flip 20" Wide | 1        |
| 1c, 2h, 3-19 | 502399      | Footrest Assembly One Piece Flip 21" Wide | 1        |
| 1c, 2i, 3-19 | 502400      | Footrest Assembly One Piece Flip 22" Wide | 1        |
| 1c, 2j, 3-19 | 114298      | Footrest Assembly One Piece Flip 23" Wide | 1        |
| 1c, 2k, 3-19 | 114299      | Footrest Assembly One Piece Flip 24" Wide | 1        |
| 1c, 2l, 3-19 | 114300      | Footrest Assembly One Piece Flip 25" Wide | 1        |
| 1c, 2m, 3-19 | 114301      | Footrest Assembly One Piece Flip 26" Wide | 1        |
| 1a           | 000226      | Footrest Platform Angle Adjustable 8"     | 1        |
| 1b           | 000227      | Footrest Platform Angle Adjustable 10"    | 1        |
| 1c           | 001319      | Footrest Platform Angle Adjustable 12"    | 1        |
| 2a           | 000357      | 14" Footrest Tube, One Piece Flip Up      | 1        |
| 2b           | 000358      | 15" Footrest Tube, One Piece Flip Up      | 1        |
| 2c           | 000359      | 16" Footrest Tube, One Piece Flip Up      | 1        |
| 2d           | 000360      | 17" Footrest Tube, One Piece Flip Up      | 1        |
| 2e           | 000361      | 18" Footrest Tube, One Piece Flip Up      | 1        |
| 2f           | 000362      | 19" Footrest Tube, One Piece Flip Up      | 1        |
| 2g           | 000363      | 20" Footrest Tube, One Piece Flip Up      | 1        |
| 2h           | 002424      | 21" Footrest Tube, One Piece Flip Up      | 1        |
| 2i           | 002425      | 22" Footrest Tube, One Piece Flip Up      | 1        |
| 2j           | 112890      | 23" Footrest Tube, One Piece Flip Up      | 1        |
| 2k           | 112891      | 24" Footrest Tube, One Piece Flip Up      | 1        |
| 2l           | 112892      | 25" Footrest Tube, One Piece Flip Up      | 1        |
| 2m           | 112893      | 26" Footrest Tube, One Piece Flip Up      | 1        |
| 3            | 100789      | Footrest Adjustable Clamp 5"              | 1        |
| 4            | 000355      | Flip Up Footrest Pivot                    | 2        |

| Item Number | Part Number | Description                               | Quantity |
|-------------|-------------|---|----------|
| 5           | 000053      | Saddle 3/4"                               | 3        |
| 6           | 100797      | M6 x 40 BHCS BLK ZC                       | 1        |
| 7           | 101455      | M6 Belleville Washer (.252 x .551 x .050) | 1        |
| 8           | 100746      | M6 Flat Washer Blk Zn                     | 1        |
| 9           | 101456      | M6 Ovl Top L/N ZC                         | 1        |
| 10          | 100669      | M5 x 10 BTN SHCS BLKZC w/Patch            | 8        |
| 11          | 000353      | Saddle - 3/4" .266 ID                     | 6        |
| 12          | 000073      | M5 x 24 Threaded Barrel                   | 4        |
| 13          | 100658      | M6 Nylock Nut Blk Zn                      | 3        |
| 14          | 100662      | M6 x 25 FSHCS BLK ZC                      | 2        |
| 15          | 000356      | Flip Up Footrest Latch                    | 1        |
| 16          | 101811      | M6 x 35 BHCS Blk Zn                       | 1        |
| 17          | 100831      | 1/2" Washer OD .040" Thick                | 2        |
| 18          | 000354      | Flip Up Ferrule                           | 1        |
| 19          | 100638      | Tube End Plug 3/4"                        | 2        |