


(Discontinued) Liberty Mag Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and Handrim kits will be listed. Handrims are not available on 12" or 16" rear wheel sizes. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed.

Item Number	Part Number	Description	Quantity
	116077	20" Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	116078	20" Mag, Full Poly Tire, Aluminum Anodized Handrim w/ Micro Grip Kit, Close Mount	1
	116079	20" Mag, Full Poly Tire, Plastic Coated Anodized Handrim Kit, Close Mount	1

Ki Mobility Mag Wheel/ Tire/ Handrim Kit Item Number Reference

Rear Tire Type	Handrim	Wheel		
		20" 5-Spoke X Core Mag	22" 5-Spoke X Core Mag	24" 5-Spoke X Core Mag
1 3/8" Street Pneumatic	Aluminum Anodized	502090	502100	502110
	Plastic Coated	502091	502101	502111
	Projection			502112
	Natural Fit Standard Omit Grip			502113
	Natural Fit Standard w/ Regular Grip			502114
	Natural Fit Standard w/ Super Grip			502115
	Natural Fit LT Omit Grip			502116
	Natural Fit LT w/ Regular Grip			502117
	Natural Fit LT w/ Super Grip			502118
1 3/8" Full Profile Poly	Aluminum Anodized	502094	502104	502127
	Plastic Coated	502095	502105	502128
	Projection			502129
	Natural Fit Standard Omit Grip			502130
	Natural Fit Standard w/ Regular Grip			502131
	Natural Fit Standard w/ Super Grip			502132
	Natural Fit LT Omit Grip			502133
	Natural Fit LT w/ Regular Grip			502134
	Natural Fit LT w/ Super Grip			502135
1 3/8" Pneumatic w/ Airless Insert	Aluminum Anodized	502096	502106	502136
	Plastic Coated	502097	502107	502137
	Projection			502138
	Natural Fit Standard Omit Grip			502139
	Natural Fit Standard w/ Regular Grip			502140
	Natural Fit Standard w/ Super Grip			502141
	Natural Fit LT Omit Grip			502142
	Natural Fit LT w/ Regular Grip			502143
	Natural Fit LT w/ Super Grip			502144

Ki Mobility Mag Wheel/ Tire/ Handrim Kit Item Number Reference

Rear Tire Type	Handrim	Wheel		
		20" 5-Spoke X Core Mag	22" 5-Spoke X Core Mag	24" 5-Spoke X Core Mag
1 3/8" Euro Trek Knobby	Aluminum Anodized			502153
	Plastic Coated			502154
	Projection			502155
	Natural Fit Standard Omit Grip			502156
	Natural Fit Standard w/ Regular Grip			502157
	Natural Fit Standard w/ Super Grip			502158
	Natural Fit LT Omit Grip			502159
	Natural Fit LT w/ Regular Grip			502160
	Natural Fit LT w/ Super Grip			502161
1" Low Profile Poly	Aluminum Anodized	502092	502102	502119
	Plastic Coated	502093	502103	502120
	Projection			
	Natural Fit Standard Omit Grip			502121
	Natural Fit Standard w/ Regular Grip			502122
	Natural Fit Standard w/ Super Grip			502123
	Natural Fit LT Omit Grip			502124
	Natural Fit LT w/ Regular Grip			502125
	Natural Fit LT w/ Super Grip			502126
1" Sentinel Solid Tire	Aluminum Anodized			
	Plastic Coated			
	Projection			
	Natural Fit Standard Omit Grip			
	Natural Fit Standard w/ Regular Grip			
	Natural Fit Standard w/ Super Grip			
	Natural Fit LT Omit Grip			
	Natural Fit LT w/ Regular Grip			
	Natural Fit LT w/ Super Grip			