

Back Tubes for Anti-Rollback System



| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| | 113299 | Backtube Straight w/ Push Handle Short with Backcane Extension Catalyst | 2 |
| | 113300 | Backtube Straight w/ Push Handle Medium with Backcane Extension, Catalyst | 2 |
| | 113301 | Backtube Straight w/ Push Handle Tall with Backcane Extension, Catalyst | 2 |
| | 113302 | Backtube 8° Bend w/ Push Handle Medium with Backcane Extension, Catalyst | 2 |
| | 113303 | Backtube 8° Bend w/ Push Handle Tall with Backcane Extension Catalyst | 2 |
| | 113304 | Backtube Straight Half Folding Back Upper w/ Push Handle Short with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Half Folding Backrest is not already present, you will need this kit along with backrest hinge latch assemblies and lower tube found on the Catalyst Half Folding Backrest page.</i> | 2 |
| | 113305 | Backtube Straight Half Folding Back Upper w/ Push Handle Tall with Backcane Extension, Catalyst | 2 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | | <i>This is just the upper backtube assembly. If the Half Folding Backrest is not already present, you will need this kit along with backrest hinge latch assemblies and lower tube found on the Catalyst Half Folding Backrest page.</i> | |
| | 113306 | Backtube 8° Bend Half Folding Back Upper w/ Push Handle Short with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Half Folding Backrest is not already present, you will need this kit along with backrest hinge latch assemblies and lower tube found on the Catalyst Half Folding Backrest page.</i> | 2 |
| | 113307 | Backtube 8° Bend Half Folding Back Upper w/ Push Handle Tall with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Half Folding Backrest is not already present, you will need this kit along with backrest hinge latch assemblies and lower tube found on the Catalyst Half Folding Backrest page.</i> | 2 |
| | 113315 | Backtube Angle Adjustable Upper w/ Push Handle Short with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Angle Adjustable Backrest is not already present, you will need this kit along with the Angle Adjustable Back Post Lower Assembly found on the Catalyst Angle Adjustable Backrest page.</i> | 2 |
| | 113308 | Backtube Angle Adjustable Upper w/ Push Handle Medium with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Angle Adjustable Backrest is not already present, you will need this kit along with the Angle Adjustable Back Post Lower Assembly found on the Catalyst Angle Adjustable Backrest page.</i> | 2 |
| | 113309 | Backtube Angle Adjustable Upper w/ Push Handle Tall with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Angle Adjustable Backrest is not already present, you will need this kit along with the Angle Adjustable Back Post Lower Assembly found on the Catalyst Angle Adjustable Backrest page.</i> | 2 |
| | 113316 | Backtube Depth Adj Height Adj Upper w/ Push Handle Short with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Height Adjustable Straight with Push Handle is not already present, you will need this kit along with retrofitting to the Catalyst Height Adjustable Straight with Push Handle (Depth Adjustable) Backrest (Kit 118429). Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i> | 2 |
| | 113310 | Backtube Depth Adj Height Adj Upper w/ Push Handle Medium with Backcane Extension, Catalyst <i>This is just the upper backtube assembly. If the Height Adjustable Straight with Push Handle is not already present, you will need this kit along with retrofitting to the Catalyst Height Adjustable Straight with Push Handle (Depth Adjustable) Backrest (Kit 118430). Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i> | 2 |
| | 113311 | Backtube Depth Adj Height Adj Upper w/ Push Handle Tall with Backcane Extension, Catalyst | 2 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | | <p><i>This is just the upper backtube assembly. If the Height Adjustable Straight with Push Handle is not already present, you will need this kit along with retrofitting to the Catalyst Height Adjustable Straight with Push Handle (Depth Adjustable) Backrest (Kit 118431). Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i></p> | |
| | 506811 | <p>Stroller Handle Back Post with Grip Assembly, Short</p> <p><i>If chair was invoiced with backrest code 50BOP032, new backcanes are not required to retrofit the Anti-Rollback System. If chair was invoiced with backrest code 50BOP010, a backrest retrofit to the new Depth Adjustable Back design is required in order to retrofit the Anti-Rollback System. See Catalyst Stroller Handle Back Post (Depth Adjustable) page. Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i></p> | 2 |
| | 506812 | <p>Stroller Handle Back Post with Grip Assembly, Tall</p> <p><i>If chair was invoiced with backrest code 50BOP032, new backcanes are not required to retrofit the Anti-Rollback System. If chair was invoiced with backrest code 50BOP010, a backrest retrofit to the new Depth Adjustable Back design is required in order to retrofit the Anti-Rollback System. See Catalyst Stroller Handle Back Post (Depth Adjustable) page. Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i></p> | 2 |
| | 112646 | <p>Stroller Handle Back Post for Removable Stroller Handle Extension, Short</p> <p><i>If chair was invoiced with backrest code 50BOP032, new backcanes are not required to retrofit the Anti-Rollback System. If chair was invoiced with backrest code 50BOP010, a backrest retrofit to the new Depth Adjustable Back design is required in order to retrofit the Anti-Rollback System. See Catalyst Stroller Handle Back Post with Removable Stroller Handle Extension (Depth Adjustable) page. Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i></p> | 2 |
| | 112647 | <p>Stroller Handle Back Post for Removable Stroller Handle Extension, Tall</p> <p><i>If chair was invoiced with backrest code 50BOP032, new backcanes are not required to retrofit the Anti-Rollback System. If chair was invoiced with backrest code 50BOP010, a backrest retrofit to the new Depth Adjustable Back design is required in order to retrofit the Anti-Rollback System. See Catalyst Stroller Handle Back Post with Removable Stroller Handle Extension (Depth Adjustable) page. Also used with the Seating Dynamics Rocker Back for full assembly please visit the Seating Dynamics Rocker Back page.</i></p> | 2 |