


## Focus Spoke Wheel / Tire / Handrim Kits



*Not all components of the Wheel, Tire and Handrim kits will be listed. These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. For Focus chairs, new wheel locks might be needed when changing wheel size or style.*

Item Number	Part Number	Description	Quantity
	118748	24" Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, 5 Rivet Nut, Standard Mount	1
	118749	24" Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, 5 Rivet Nut, Close Mount	1
	118752	24" Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Standard Mount	1
	118753	24" Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Close Mount	1
	118756	24" Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, 5 Rivet Nut, Standard Mount	1
	118757	24" Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, 5 Rivet Nut, Close Mount	1
	118760	24" Spoke, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118761	24" Spoke, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount	1
	118978	24" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118979	24" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118980	24" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip, Standard Mount	1
	118981	24" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip, Close Mount	1
	118800	24" Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, 5 Rivet Nut, Standard Mount	1
	118801	24" Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, 5 Rivet Nut, Close Mount	1
	118804	24" Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Standard Mount	1
	118805	24" Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Close Mount	1
	118808	24" Spoke, Full Poly Tire, Plastic Coated Handrim Kit, 5 Rivet Nut, Standard Mount	1
	118809	24" Spoke, Full Poly Tire, Plastic Coated Handrim Kit, 5 Rivet Nut, Close Mount	1
	118812	24" Spoke, Full Poly Tire, Projection Handrim Kit, Standard Mount	1
	118813	24" Spoke, Full Poly Tire, Projection Handrim Kit, Close Mount	1